

Waterloo Healthcare



Contract Holder
Contract V797P-4115B

Anesthesia Carts



Page 16

Crash Carts



Page 6

Isolation Carts



Page 18

Bedside/Slim Carts



Page 20

Treatment/Procedure Carts



Page 22

Medication Carts



Page 43

Toll Free 1-800-833-4419 Direct 602-414-3691
www.waterloohealthcare.com



My Account Number is: _____

My Local Sales Representative:

TABLE OF CONTENTS



NEW PRODUCTS

- ❖ Cassette Cart4
- ❖ Hanging Organizers4
- ❖ Proximity Automated Locking Systems5

CART HIGHLIGHTS

- ❖ Lightweight Aluminum Carts6-7
- ❖ Traditional Steel Carts8-9

EMERGENCY CARTS AND PACKAGES

- ❖ Instant Access Carts (ER2000).....10
- ❖ Emergency Carts & Packages11-14
- ❖ Pediatric Carts & Packages.....15

ANESTHESIA CARTS

- ❖ Lightweight Aluminum Anesthesia Carts16
- ❖ Traditional Steel Anesthesia Carts17

ISOLATION / INFECTION CONTROL CARTS

- ❖ Lightweight Aluminum Isolation Carts18
- ❖ Traditional Steel Isolation Carts19

HANGING ORGANIZERS18

SLIM / BEDSIDE CARTS20-21

TREATMENT AND PROCEDURE CARTS

- ❖ Airway22
- ❖ Split Drawer Cart.....22
- ❖ Cast.....23
- ❖ IV Therapy.....23
- ❖ Simulated Wood24
- ❖ Mobile Workstation24
- ❖ Phlebotomy.....25
- ❖ Economy Series.....26-27

CUSTOMIZE YOUR OWN CART28-31



TABLE OF CONTENTS



continued...

ACCESSORIES AND DIVIDERS

- ❖ Unicart Accessory Options32-36
- ❖ Additional Locking and Security37
- ❖ Trays and Divider Systems38-40
- ❖ Accessory Packages.....40
- ❖ Portable Storage41

NARCOTICS CABINETS42

MEDICATION CART HIGHLIGHTS

- ❖ Simulated Wood Cart.....44
- ❖ MEDTRX Series Medication Carts46
- ❖ Dual-Sided Medication Carts47

PATIENT CASSETTE CARTS

- ❖ Half-Sized Cassette Carts.....48
- ❖ Full-Sized Cassette Carts48
- ❖ Dual-Sided Cassette Carts49

UNIT DOSE BOX CARTS50

PILL CARD CARTS50-51

BINS, CASSETTES, ACCESSORIES

- ❖ Exchange Cassettes and Bins52-53
- ❖ Pill Card Storage Bins and Dividers54
- ❖ Medication Cart Accessories.....55-56

CUSTOM MEDICATION CART CONFIGURATION57

TERMS AND CONDITIONS.....58

INDEX59-61

COLOR OPTIONS62

SPECIFICATIONS.....63

New Products

Hanging Organizers

Constructed with durable & lightweight aluminum.

HANGORG-1



Dimensions:
20"W x 26.75"H x 4.5"D

Holds:
4 Glove Boxes
2 Face Mask box
Flat Pack Gowns

Mounting:

- Permanent mounting with screws (not included).
- Hanging from doorway with optional accessory (ISOHANG-1).

ISOORG-1



Dimensions:
16.5"W x 25.25"H x 4.5"D

Holds:
3 Glove Boxes
1 Face Mask box
Flat Pack Gowns

Mounting:

- Permanent mounting with screws (not included).
- Hanging from doorway with optional accessory (ISOHANG-1).

ISOORG-2



Dimensions:
22.5"W x 20.5"H x 4.75"D

Holds:
2 Glove Boxes
1 Face Mask Box
1 Face Shield Box
Flat Pack Gowns

Mounting:

- Permanent mounting with screws (not included).
- Hanging from doorway with optional accessory (ISOHANG-1).


Cassette Cart


Designed to store essential equipment, materials and medications while keeping you organized. The cassette locks separately from the cart allowing access only to authorized individuals.



UTGKU-C3336-LTG

Features & Benefits

- Antimicrobial properties built-in to protect portions of the product.  **Antimicrobial Surfaces**
- Centralized locking with 3 available options
- Drawers with ball bearing slides with quick release
- Available in 16 drawer colors (see below)
- Dual push handles
- ABS plastic molded top.
- 5" Premium easy roll swivel casters (2 brake, 1 tracking, 1 full swivel).
- Stabilizer base with high impact plastic bumper to add stability
- Latex Free
- Bins sold separately, choice of 3", 4", 6" wide bins
- Available in three sizes:
Short (22" drawer height)
 29" W x 24.5" D x 37.5" H
Tall (28" drawer height)
 29" W x 24.5" D x 43.5" H
X-Tall (31" drawer height)
 29" W x 24.5" D x 46.5" H



Larger View – See Page 62

New Products

Proximity Automated Locking Systems

Proximity Locking System Features & Benefits

- Available on our Lightweight Aluminum and Slim Med Jr. Unicarts.
- Allows user to use a proximity card, or a combination to gain access to the cart.
- Keypad can be disabled easily so only proximity cards are accepted.
- Uses industry-standard HID proximity system.
- Can be used with hospital's existing systems. We also offer 2 different types of proximity cards for purchase (standard credit card size and miniature size).
- Holds up to 250 user or supervisor codes and can track the last 1,500 successful and failed access attempts.
- Manual key override for emergency access.
- Auto Lock & Unlock. Ability to relock after designated amount of time (set by user). Can also be disabled if preferred.
- LED indicators show the battery status (locked or unlocked).
- Windows PC optional software package enables you to manage single or multiple carts, user codes and auditing capabilities.

The new Proximity locking systems allows users to access essential supplies in seconds. Simply swipe your badge in front of the reader and the fully automated system unlocks instantly.



UTGPA-333369-ELB


5" Premium easy roll casters with injection molded abs stabilizer base

- (2) Braking
- (1) Tracking
- (1) Total Swivel


Use the manual key override for emergency access.

Keypad holds up to 250 user or supervisor codes.

Simply wave your badge to unlock instantly.



Larger View – See Page 62



Lightweight Aluminum Cart Highlights

Standard Features & Benefits

Our ELITE Aluminum Cart offers easier maneuverability and safety when time is of the essence. Lightweight aluminum construction is 30% lighter than standard steel carts.



Antimicrobial Surfaces

- The patented U-Divide™ System included in every drawer!
- Antimicrobial properties built-in to protect portions of the product.
- New dual pull-out shelves provide additional work surface on both sides of the cart.
- 4 vertical aluminum corner mounting tracks to mount accessories without the need of predrilled holes.
- 2 side mounting panels to assist in placement of accessories.
- Dual push handles.
- Exterior comes standard in three colors: Red, Blue or Light Gray.
- ABS cutaway plastic top minimizes wrist strain when writing or working on the cart's surface.
- Centralized locking with 5 available options.
- Ability to configure drawers to fit individual needs at no additional charge.
- 5" Premium easy roll casters. (1) full swivel, (2) swivel w/ brake, (1) swivel w/ tracking.
- Stabilizer Base with high impact, non-marking plastic bumper to add stability and help keep cart from tipping.
- Drawers with ball bearing slides were designed to push in and lock even after the cart has been secured.



**Up to 40%
Lighter Weight**
than comparable steel carts.

**Drawer Dividers
INCLUDED!**
\$400
average savings!

Flexible Security Options

Let us increase the security of your cart! Gatelocks can be easily added to any cart. Some of the many uses include:

- Storing Narcotics
- Visible Inventory Management
- Securing Liquids

See page 37 for more details.

Locking Options



Gatelock




Pushbutton Locking



Key Locking



Breakaway Locking



Proximity Locking



Electronic Locking

Lightweight Aluminum Cart Features

How to customize your Waterloo Aluminum Cart

1. Choose Frame Size
2. Choose Locking Mechanism - See Page 6 to review locking options
3. Choose Drawer Front Colors - See Page 62 for color options



(To help you choose the right cart for you, answer the following questions: What size of cart do you need & what size drawers do you want?)

Waterloo offers 4 different size frames for our Lightweight Aluminum cart:

Short Aluminum Series

USGLA-399-YEL



USGLA-3369-TLG



USGLA-33366-RED



3-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 36"H
 (includes stabilizer system and casters)
Drawer Sizes: 22"W x 16.5"D
 21" Configurable drawer space
 (1) 3"H, (2) 9"H

4-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 36"H
 (includes stabilizer system and casters)
Drawer Sizes: 22"W x 16.5"D
 21" Configurable drawer space
 (2) 3"H, (1) 6"H, (1) 9"H

5-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 36"H
 (includes stabilizer system and casters)
Drawer Sizes: 22"W x 16.5"D
 21" Configurable drawer space
 (3) 3"H, (2) 6"H

Medium Aluminum Series

UMGLA-33612-YEL



UMGLA-6612-LTB



UMGLA-33369-RED



4-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 39"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 24" Configurable drawer space
 (2) 3"H, (1) 6"H, (1) 12"H

3-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 39"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 24" Configurable drawer space
 (2) 6"H, (1) 12"H

5-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 39"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 24" Configurable drawer space
 (3) 3"H, (1) 6"H, (1) 9"H

Tall Aluminum Series

UTGKA-3699-LTG



UTGKA-33669-VIL



UTGLA-333369-RED



4-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 42"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 27" Configurable drawer space
 (1) 3"H, (1) 6"H, (2) 9"H

5-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 42"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 27" Configurable drawer space
 (2) 3"H, (2) 6"H, (1) 9"H

6-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 42"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 27" Configurable drawer space
 (4) 3"H, (1) 6"H, (1) 9"H

Extra-Tall Aluminum Series

UXGKA-36669-LTB



UXGKA-333669-NTL



UXGKA-3333369-FWG



5-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 45"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 30" Configurable drawer space
 (1) 3"H (3) 6"H, (1) 9"H

6-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 45"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 30" Configurable drawer space
 (3) 3"H, (2) 6"H, (1) 9"H

7-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 45"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 30" Configurable drawer space
 (5) 3"H, (1) 6"H, (1) 9"H

Aluminum Carts

Traditional Steel Cart Highlights

Standard Features:

- Antimicrobial properties built-in to protect portions of the product.
- Ball bearing slides with quick release.
- Dual push handles.
- ABS plastic molded top.
- 5" Premium easy roll casters.
- Stabilizer Base with high impact, non-marking plastic bumper to add stability and help keep cart from tipping.
- Pre-threaded holes on sides & back to easily mount accessories.
- Latex Free.



Antimicrobial Surfaces

Breakaway Locking Cart

With our centralized locking you will never have to reach down to unlock the cart again!

Just a simple flick of your wrist and the lever will break the seal and allow access to every drawer.



USRLU-4369-RED

Key Locking Cart



UTGKU-4699-YEL

5" Premium easy roll casters

- (1) Swivel
- (2) Swivel w/Brake
- (1) Tracking Caster

3 sizes available to fit your needs

Short - 29"W x 24.5"D x 37.5"H

Tall - 29"W x 24.5"D x 43.5"H

X-Tall - 29"W x 24.5"D x 46.5"H

Push Button Locking Cart

The push button locking cart comes standard with lever lock and override key.

Stabilizer Base with high impact plastic bumper to add stability by lowering the center of gravity on the cart, keeping it from tipping.



UTGSU-433369-DKB

Traditional Steel Cart Features

How to customize your Waterloo Steel Unicart

1. Choose Frame Size
2. Choose Locking Mechanism - See Page 8 to review locking options
3. Choose Drawer Front Colors - See below for color options



(To help you choose the right cart for you, answer the following questions: What size of cart do you need & what size drawers do you want?)

Waterloo offers 3 different size frames for our traditional steel cart:

Short Uni-Cart Series



3-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 37.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 22" Configurable drawer space
 (1) 4"H, (2) 9"H

4-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 37.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 22" Configurable drawer space
 (1) 4"H, (1) 3"H, (1) 6"H, (1) 9"H

5-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 37.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 22" Configurable drawer space
 (1) 4"H, (2) 3"H, (2) 6"H

Tall Uni-Cart Series



4-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 43.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 28" Configurable drawer space
 (1) 4"H, (1) 6"H, (2) 9"H

5-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 43.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 28" Configurable drawer space
 (1) 4"H, (1) 3"H, (2) 6"H, (1) 9"H

6-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 43.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 28" Configurable drawer space
 (1) 4"H, (3) 3"H, (1) 6"H, (1) 9"H

Extra-Tall Uni-Cart Series




5-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 46.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 31" Configurable drawer space
 (1) 4"H, (3) 6"H, (1) 9"H

6-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 46.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 31" Configurable drawer space
 (1) 4"H, (2) 3"H, (2) 6"H, (1) 9"H

7-Drawer Cart Dimensions
Overall: 29"W x 24.5"D x 46.5"H
 (includes stabilizer system and casters)
Drawer Sizes 22"W x 16.5"D
 31" Configurable drawer space
 (1) 4"H, (4) 3"H, (1) 6"H, (1) 9"H

Larger View – See Page 62



Steel Carts

ER2000™ Instant Access Emergency Cart

This emergency cart is unlike any other on the market. For your emergency department, standardize on the ER2000™ series. A state-of-the-art design that is attractive yet functional. Simple and easy to use, press the bright yellow “open” button located at the top center. The doors spring open and can lock into place allowing full access to the contents needed in any emergency situation.



Features & Benefits

- Antimicrobial properties built-in to protect portions of the product.
- Instant access – one touch breakaway seal opens for quick response.
- Door latches keep doors in open stationary position.
- Full extension see-thru drawer trays
- Recessed O2 tank brackets protect valves and gauges.
- Removable transportable trays for ease of restocking with sealable bag system.
- Quiet, smooth, ball bearing slides for years of reliable drawer operation.
- Adjustable door storage for more flexibility.
- Additional pull out work surface.
- Tough molded stabilizer frame bumper allows all drawers to be fully extended without risk of tipping cart.
- 5” Premium easy roll casters (two brake, one tracking, one swivel)

Included Accessories

- Cardiac board (CBB-ER)
- IV Pole (IV-3)
- (2) Ambu bag hooks (ABH-1)
- Plastic breakaway seals (PS-100)
- Oxygen tank brackets (O-2)
- Airway, Breathing, Circulation Labels (ER-LBL)

Instant Access Cart

Dimensions (includes casters, bumper and top)


Overall.....26”W x 33”D x 45”H
 Open43”W x 33”D x 45”H

Drawer trays

(1-4)22.25”W x 16.5”D x 2.75”H
 (5).....22.25”W x 16.5”D x 5.5”H

Door bins

(4 each)7.5”W x 2.9”D x 2.8”H
 (2 each)7.2”W x 5.6”D x 2.6”H



Additional Surface Area

Sliding Work Surface

extends 11.25” from either side for additional work area.

Complete Emergency Cart Packages

Lightweight Aluminum Emergency Cart Package



EMER-PKGA-DLX

Deluxe Emergency/Crash accessory package includes:

- Cardiac board (CBB-7)
- IV pole and bracket (IV-2)
- Oxygen tank brackets (OTB-4)
- Bag of (100) plastic seals (PS-100)
- Railing and supports (R-4)
- Base rail system (BRS-3)
- Top shelf (SHLF-3),
- (2) 3" Divisible trays with dividers (WDC-2025)
- (2) 5" Divisible trays with dividers (WDC-2050)
- 3" drawer tray with lid (EMER-BOX)
- Dividers for EMER-BOX (EMER-DV)
- Hospital grade outlets strip (OS-1)
- Suction shelf (SS-M)

Traditional Steel Emergency Cart Package

EMER-PKG-DLX

Deluxe Emergency/Crash accessory package includes:

- Cardiac board (CBB-6)
- IV pole and bracket (IV-2)
- Oxygen tank brackets (OTB-4)
- Bag of (100) plastic seals (PS-100)
- Collapsible side shelf (SSG-2C)
- (2) utility hooks (UH-2)
- Base rail system (BRS-2)
- Top shelf (SHLF-2)
- Railing and supports (R-3)
- (2) 3" Divisible trays with dividers (WDC-2025)
- (2) 5" Divisible trays with dividers (WDC-2050)
- 3" drawer tray with lid (EMER-BOX)
- Dividers for EMER-BOX (EMER-DV)
- Hospital grade outlet strip (OS-1)
- Suction shelf (SS-M)



Standard Emergency Cart Packages

Lightweight Aluminum Emergency Cart Package

EMER-PKGA

Standard Emergency / Crash accessory package for lightweight aluminum Unicarts includes:

(shown at right, cart not included)

- Cardiac board (CBB-7)
- IV pole and bracket (IV-2)
- Oxygen tank brackets (OTB-4)
- Bag of (100) plastic seals (PS-100)
- (2) Utility hooks (UH-2)
- Railing and supports (R-4)
- (2) 3" Divisible trays with dividers (WDC-2025)
- (2) 5" Divisible trays with dividers (WDC-2050)



Antimicrobial Surfaces

Dual slide-out shelves.



Drawer Dividers INCLUDED!

\$400

average savings!

USRLA-3369-RED



4-Drawer Cart Dimensions

Overall: 29"W x 24.5"D x 36"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
21" Configurable drawer space
(2) 3"H, (1) 6"H, (1) 9"H

Standard Cart Pkgs



UTRLU-43336P-RED

Can be configured as 9" drawer also.

5-Drawer w/Panel Cart Dimensions

Overall: 29"W x 24.5"D x 43.5"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
31" Configurable drawer space
(1) 4"H, (3) 3"H, (1) 6"H, (1) 9" Panel

Traditional Steel Emergency Cart Package

EMER-PKGM

Traditional Steel Unicart with Standard Emergency / Crash accessory package includes: (shown at left, cart not included)

- Cardiac board (CBB-6)
- IV pole and bracket (IV-2)
- Oxygen tank brackets (OTB-4)
- Bag of (100) plastic seals (PS-100)
- Collapsible side shelf (SSG-2C)
- (2) utility hooks (UH-2)
- Railing and supports (R-3)
- 4 x clear plastic trays with dividers (CP-DIVTRAY)
- 2" cut and clip divider system (DIV-CC2)

Single Colored Emergency Carts

Lightweight Aluminum



Antimicrobial Surfaces



6-Drawer Cart Dimensions UXRLA-333669-RED

Overall: 29"W x 24.5"D x 45"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
27" Configurable drawer space
(3) 3"H, (2) 6"H, (1) 9"H

5-Drawer Cart Dimensions UTRLA-33669-RED

Overall: 29"W x 24.5"D x 42"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
27" Configurable drawer space
(2) 3"H, (2) 6"H, (1) 9"H

4-Drawer Cart Dimensions UMRLA-3669-RED

Overall: 29"W x 24.5"D x 39"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
21" Configurable drawer space
(1) 3"H, (2) 6"H, (1) 9"H

4-Drawer Cart Dimensions USRLA-3369-RED

Overall: 29"W x 24.5"D x 36"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
21" Configurable drawer space
(2) 3"H, (1) 6"H, (1) 9"H

Single Colored Carts

Traditional Steel



6-Drawer Cart Dimensions UXRLU-433669-RED

Overall: 29"W x 24.5"D x 43.5"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
28" Configurable drawer space
(1) 4"H, (3) 3"H, (1) 6"H, (1) 9"H Panel



5-Drawer Cart Dimensions UTRLU-43669-RED

Overall: 29"W x 24.5"D x 43.5"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
28" Configurable drawer space
(1) 4"H, (1) 3"H, (2) 6"H, (1) 9"H



4-Drawer Cart Dimensions USRLU-4369-RED

Overall: 29"W x 24.5"D x 37.5"H
(includes stabilizer system and casters)
Drawer sizes: 22"W x 16.5"D
22" Configurable drawer space
(1) 4"H, (1) 3"H, (1) 6"H, (1) 9"H

Larger View – See Page 62



GRA-Charcoal



DKB-Dark Blue



ELB-Electric Blue



FWG-Fairway Green



LTB-Light Blue



LTG-Light Gray



NTL-Natural



ORG-Orange



PNK-Pink



RED-Red




SPG-Spring Green



TLG-Teal



VIL-Violet



WHT-White



WNE-Wine



YEL-Yellow

Multicolored Emergency Carts



Multicolored Carts

Aluminum



6-Drawer Cart Dimensions
UTGLA-333369-ELB
 Overall: 29"W x 24.5"D x 42"H
 (includes stabilizer system and casters)
 Drawer sizes: 22"W x 16.5"D
 27" Configurable drawer space
 (4) 3"H, (1) 6"H, (1) 9"H

Steel



6-Drawer Cart Dimensions
UTGLU-433369-RED
 Overall: 29"W x 24.5"D x 42"H
 (includes stabilizer system and casters)
 Drawer sizes: 22"W x 16.5"D
 28" Configurable drawer space
 (1) 4"H, (3) 3"H, (1) 6"H, (1) 9"H

Space Saver



4-Drawer Cart Dimensions
JSGLU-3369-ELB
 Overall: 23"W x 24.5"D x 36"H
 (includes stabilizer system and casters)
 Drawer sizes: 16.25"W x 16.5"D
 21" Configurable drawer space
 (2) 3"H, (1) 6"H, (1) 9"H




6-Drawer Cart Dimensions
UXGLA-333669-RED
 Overall: 29"W x 24.5"D x 45"H
 (includes stabilizer system and casters)
 Drawer sizes: 22"W x 16.5"D
 30" Configurable drawer space
 (3) 3"H, (2) 6"H, (1) 9"H



5-Drawer Cart with Panel Dimensions
UTGLU-43336P-ELB
 Overall: 29"W x 24.5"D x 42"H
 (includes stabilizer system and casters)
 Drawer sizes: 22"W x 16.5"D
 28" Configurable drawer space
 (1) 4"H, (3) 3"H, (1) 6"H, (1) 9" Panel



6-Drawer Cart Dimensions
JTGLU-333369-RED
 Overall: 23"W x 24.5"D x 42"H
 (includes stabilizer system and casters)
 Drawer sizes: 16.25"W x 16.5"D
 27" Configurable drawer space
 (4) 3"H, (1) 6"H, (1) 9"H



Larger View – See Page 62

Pediatric Emergency Carts



Lightweight Aluminum Pediatric Carts



This pediatric cart has the same features as our Tall Elite Aluminum cart and is designed with nine color-coded drawers for your pediatric needs. Many cheerful colors to choose from to brighten any pediatric department.



Dual slide-out shelves.



Drawer Dividers INCLUDED!

\$400

average savings!

PEDS-PKGA package includes:

(shown at left, cart not included)

- Base rail system (BRS-3)
- Top Shelf (SHLF-3)
- Tilt Bins (TB-5)
- Waste Container (BSK-2)
- Railing and Supports (R-4)
- Chart Holder (CHT-1)
- IV Pole (IV-2)
- Plastic Breakaway Seals (PS-100)
- Three-inch gate lock bar (GLB-3)

UTGLA-9PEDS

Also available in 10 drawer models.

Contact customer service for product details.

PEDS-PKG package:

(shown at right, cart not included)

- Base rail system (BRS-2)
- Waste receptacle (BSK-1)
- IV pole (IV-2)
- Top shelf (SHLF-2)
- Side rail (2 each SR-1)
- Chart holder (CHT-1) (not shown)
- Collapsible side shelf (SSG-2C)
- 5 compartment tilt bins (TB-5)
- Railing and supports (R-3)

See page 32-39 for more accessories.

Traditional Steel Pediatric Carts


UTGLU-9PEDS

shown with optional PEDS-PKG



Add a gatelock to any of our Uni-carts. The locks can be customized to lock one or more drawers (see page 37).

Larger View – See Page 62



Lightweight Aluminum Anesthesia Carts

Offering accessories and flexibility for a demanding department.

ANES-PKGA-DLX


(shown at right, cart not included)

Deluxe Anesthesia "Elite" cart accessory package includes:

- Waste receptacle (BSK-1)
- Base rail system (BRS-3)
- 4 Compartment tilt bins (TB-4)
- 6 Compartment tilt bins (TB-6)
- Tape dispenser (TD-1)
- Railing and supports (R-4)
- Wire storage basket (WSB-1)
- Catheter holder (CH-3)
- Utility receptacle (UR-2)
- Utility receptacle (UR-3)
- (4) 3" Divisible trays with dividers (WDC-2025)
- (2) 5" Divisible trays with dividers (WDC-2050)
- Modular drawer tray with ampule holder (DIVTRAY-AMP)
- Clear organizer (ORG-1S)



Dual slide-out shelves.



Drawer Dividers INCLUDED!

\$400

average savings!



UTGEA-333369-DKB
w/deluxe package



UTGKA-333369-FWG w/standard package



Antimicrobial Surfaces

ANES-PKGA

(shown at left, cart not included)

Standard Anesthesia "Elite" cart accessory package includes:

- Waste receptacle (BSK-1)
- Base rail system (BRS-3)
- 6 Compartment tilt bins (TB-6)
- Tape dispenser (TD-1)
- Clear organizer (ORG-1S)
- Shown with optional wrap around railing (R-4)

Looking for customizable drawer trays? Turn to pages 38-39!



DIVTRAY-ANES

Our most popular divider set. Includes pre-configured dividers and tray. 21.5"W x 15.25"D x 2.3"H



DIVTRAY-AMP

Ampule holder divider system with removable tray. Fits in all full-size drawers. 21.5"W x 15.25"D x 2.3"H

Traditional Steel Anesthesia Carts

Offering accessories and flexibility for a demanding department.

ANES-PKG-DLX

(shown at right, cart not included)

Deluxe Anesthesia Traditional Steel cart accessory package includes:

- Waste receptacle (BSK-1)
- Base rail system (BRS-2)
- 4 Compartment tilt bins (TB-4)
- 6 Compartment tilt bins (TB-6)
- Tape dispenser (TD-1)
- Railing and supports (R-3)
- Wire storage basket (WSB-1)
- Catheter holder (CH-3)
- Utility receptacle (UR-2)
- Utility receptacle (UR-3)
- (4) 3" Divisible trays with dividers (WDC-2025)
- (2) 5" Divisible trays with dividers (WDC-2050)
- Modular drawer tray with ampule holder (DIVTRAY-AMP)
- Clear organizer (ORG-1S)
- (2) Side rails to mount accessories (SR-1)



Antimicrobial Surfaces



UTGSU-433369-FWG

ANES-PKG

(shown at left, cart not included)

Standard Anesthesia Traditional Steel cart accessory package includes:

- Waste receptacle (BSK-1)
- Base rail system (BRS-2)
- 6 Compartment tilt bins (TB-6)
- Tape dispenser (TD-1)
- Small clear organizer (ORG-1S)
- (2) Side rails to mount accessories (SR-1)

Anesthesia Carts



UTGKU-433369-DKB



WMT-3

Full drawer tray with three rows of ampule holders. Total of 33 dividers, four securing clips and removable tray. Dimensions: 21.5"W x 15.25"D x 2.3"H



DIVTRAY-STD

Standard divider system used to store larger items. Includes dividers and removable tray. Fits in all full-size drawers. 21.5"W x 15.25"D x 2.3"H

Lightweight Aluminum Carts & Organizers

Infectious disease programs are critical issues facing all hospitals. Waterloo Isolation Systems address those preventative measures for your facility.



Hanging Organizers

Constructed with durable & lightweight aluminum.

HANGORG-1



Dimensions:
20"W x 26.75"H x 4.5"D

Holds:
4 Glove Boxes
2 Face Mask box
Flat Pack Gowns

Mounting:

- Permanent mounting with screws (not included).
- Hanging from doorway with optional accessory (ISOHANG-1).

ISOORG-1



Dimensions:
16.5"W x 25.25"H x 4.5"D

Holds:
3 Glove Boxes
1 Face Mask box
Flat Pack Gowns

Mounting:

- Permanent mounting with screws (not included).
- Hanging from doorway with optional accessory (ISOHANG-1).

ISOORG-2



Dimensions:
22.5"W x 20.5"H x 4.75"D

Holds:
2 Glove Boxes
1 Face Mask Box
1 Face Shield Box
Flat Pack Gowns

Mounting:

- Permanent mounting with screws (not included).
- Hanging from doorway with optional accessory (ISOHANG-1).


INF-PKGA
(shown on right, cart not included)

Isolation/Infection accessory package for Elite Aluminum includes:

- Coat rack (CR-1)
- 2 each Utility hooks (UH-2)
- Railing and supports (R-4)

UTGKA-6669-YEL
4-drawer Elite Aluminum isolation cart
Drawers: (3) 6"H, (1) 9"H
Accessory package not included.

For Elite cart features and dimensions see page 6 and 7.



UTGKA-999-YEL
3-drawer Elite Aluminum isolation cart
Drawers: (3) 9"H
Available in five different locking options. See page 6.

UXGKA-3999-YEL
4-drawer Elite Aluminum isolation cart
Drawers: (1) 3"H, (3) 9"H
For Elite cart features and dimensions see page 6 and 7.



Larger View – See Page 62

Isolation Carts



Traditional Steel Series

Antimicrobial properties are built in to protect portions of ALL Waterloo Healthcare carts.



USGKU-499-YEL

3-drawer Traditional Steel Isolation cart. Drawers: (1) 4”H, (2) 9”H

For Steel cart features and dimensions see page 8 and 9.



UTGKU-4699-YEL

4-drawer traditional steel isolation cart. Drawers: (1) 4”H, (1) 6”H, (2) 9”H

Shown with optional isolation/infection package (INF-PKGM). Available in three different locking options.



UXGKU-4999-YEL

4-drawer Traditional Steel Isolation cart. Drawers: (1) 4”H, (3) 9”H

For Steel cart features and dimensions see page 8 and 9.



JSGKU-399-YEL

3-drawer Slim Med Junior Isolation Cart. Drawers: (1) 3”H, (2) 9”H

Isolation/Infection accessory package for traditional steel carts includes:

- Coat rack (CR-1)
- 2 each Utility hooks (UH-2)
- Railing and supports (R-3)

Slim Med Junior Series

For Med Junior Infection Package use part number INF-PKGJ. See page 20 for Med Jr features and benefits.

















Our Slim cart is also compatible with most accessories that fit the sides of our full-sized carts. Some drilling may be required. See pages 32-39 for options.



JTGLU-3699-YEL

4-drawer Slim Med Junior Isolation Cart. Drawers: (1) 3”H, (1) 6”H (2) 9”H

Larger View – See Page 62

| | | | | | | | |
|--|---|---|---|--|---|---|--|
|  GRA-Charcoal |  DKB-Dark Blue |  ELB-Electric Blue |  FWG-Fairway Green |  LTB-Light Blue |  LTG-Light Gray |  NTL-Natural |  ORG-Orange |
|  PNK-Pink |  RED-Red |  SPG-Spring Green |  TLG-Teal |  VIL-Violet |  WHT-White |  WNE-Wine |  YEL-Yellow |

Isolation Carts


Lightweight Aluminum Slim / Bedside Carts

Standard Features & Benefits

Our ELITE Lightweight Slim Cart offers easier maneuverability and safety when time is of the essence. Lightweight aluminum construction is 30% lighter than standard steel carts.

- The patented U-Divide™ System included in every drawer!
- Antimicrobial properties built-in to protect portions of the product.
- New dual pull-out shelves provide additional work surface on both sides of the cart.
- 4 vertical aluminum corner mounting tracks to mount accessories without the need of predrilled holes.
- 2 side mounting panels to assist in placement of accessories.
- Dual push handles.
- ABS plastic top minimizes wrist strain when writing or working on the cart's surface.
- Centralized locking with 5 available options.
- Ability to configure drawers to fit individual needs at no additional charge.
- 5" Premium easy roll casters. (1) full swivel, (2) swivel w/ brake, (1) swivel w/ tracking.
- Stabilizer Base with high impact, non-marking plastic bumper to add stability and help keep cart from tipping.
- Drawers with ball bearing slides were designed to push in and lock even after the cart has been secured.

Dual slide-out shelves.



Drawer Dividers INCLUDED!

\$400

average savings!

JTGLA-333369-ELB

Stylish drawer pulls.



JTGLA-333369-YEL



Antimicrobial Surfaces

Locking Options



Key Lock



Lever/Breakaway Lock



Pushbutton Lock



Electronic Lock



Proximity Lock



GRA-Charcoal



DKB-Dark Blue




ELB-Electric Blue



FWG-Fairway Green



LTB-Light Blue



LTG-Light Gray



NTL-Natural



ORG-Orange



PNK-Pink



RED-Red



SPG-Spring Green



TLG-Teal



VIL-Violet



WHT-White



WNE-Wine



YEL-Yellow

Larger View - See Page 62

Lightweight Aluminum Slim Med Jr Carts

JSGEA-399-DKB

JMGEA-3669-SPG

JTGEA-3699-RED

JXGEA-333669-VIL



**3-Drawer 21'' Steel
Med Jr. Emergency Cart**

JSGEA-399-DKB Light gray exterior,
Dark Blue drawers
Dimensions: 23"W x 24.5"D x 36"H
Drawer sizes: (1) 3"H, (2) 9"H

**4-Drawer 24'' Aluminum
Med Jr. Emergency Cart**

JMGEA-3669-SPG Light gray exterior,
Spring Green drawers
Dimensions: 23"W x 24.5"D x 39"H
Drawer sizes: (1) 3"H, (2) 6"H, (1) 9"H

**4-Drawer 27'' Aluminum
Med Jr. Emergency Cart**

JTGEA-3699-RED Light gray exterior,
Red drawers
Dimensions: 23"W x 24.5"D x 42"H
Drawer sizes: (1) 3"H, (1) 6"H, (2) 9"H

**6-Drawer 30'' Aluminum
Med Jr. Emergency Cart**

JXGEA-333669-VIL Light gray exterior,
Violet drawers
Dimensions: 23"W x 24.5"D x 45"H
Drawer sizes: (3) 3"H, (2) 6"H, (1) 9"H



**Antimicrobial
Surfaces**

Steel Slim Med Jr Carts

JSGKU-399-YEL

JTGKU-3999-LTG

JXGSU-333669-DKB



**3-Drawer 21'' Steel
Med Jr. Emergency Cart**

JSGKU-399-YEL Light gray exterior,
Yellow drawers
Dimensions: 23"W x 24.5"D x 36"H
Drawer sizes: (1) 3"H, (2) 9"H

**4-Drawer 27'' Steel
Med Jr. Emergency Cart**

JTGKU-4999-LTG Light gray exterior,
Light Gray drawers
Dimensions: 23"W x 24.5"D x 42"H
Drawer sizes: (1) 3"H, (3) 9"H

**6-Drawer 30'' Steel
Med Jr. Emergency Cart**

JXGSU-333669-DKB Light gray exterior,
Dark Blue drawers
Dimensions: 23"W x 24.5"D x 45"H
Drawer sizes: (3) 3"H, (2) 6"H, (1) 9"H

**Slim
Med
Junior**

Treatment and Procedure



Airway

UTGLA-333369-LTG 6-Drawer Lightweight Aluminum Unicart with Light Gray drawers

Dimensions

(includes casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 42"H

Drawers: 22"W x 16.5"D

Drawer Height: (4) 3"H, (1) 6"H, (1) 9"H

Weight: 74 lbs

This Difficult Airway cart is shown with the following optional accessories:

- (1) OS-1 – Six outlet hospital grade outlet strip
- (1) SH-2 – Two scope, side mounted scope cabinet
- (1) SHLF-110 – Two-tier shelving unit
- (1) BSK-2 – Large waste basket



Split-Drawer



UTGKU-9SR-GRA – 9-Drawer, Traditional Steel Unicart with Charcoal Gray drawers

Dimensions

(include casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 43.5"H

Drawers:


(1-2).....10.63"W x 16.5"D x 3.69"H

(3-4).....10.63"W x 16.5"D x 2.46"H

(5-7).....22"W x 16.5"D x 2.46"H

(8-9).....22"W x 16.5"D x 5.85"H

Weight: 125 lbs



Larger View – See Page 62

Treatment and Procedure



Cast

UTGLA-3333366-BLG – 7-Drawer, Lightweight Aluminum Unicart with Blue Gray drawers

Dimensions

(include casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 42"H

Drawers: 22"W x 16.5"D

Drawer Height: (5) 3"H, (2) 6"H

Weight: 79 lbs

This Cast cart is shown with the following optional accessories:

- (1) CAST-1 – Cast bucket holder and towel bar
- (1) BSK-2 – Large waste basket



IV Therapy

UTGLU-43669-PNK – 5-Drawer, Traditional Steel Unicart with Pink drawers

Dimensions

(include casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 43.5"H

Drawers: 22"W x 16.5"D

Drawer Height: (1) 4"H, (1) 3"H, (2) 6"H, (1) 9"H


Weight: 105 lbs

This IV Therapy cart is shown with the following optional accessories:

- (1) IV-2 – IV pole and bracket
- (1) SSG-2C – Collapsible side shelf
- (1) BH-2 – 2 box glove box holder
- (1) BSK-2 – Large waste basket



Larger View – See Page 62



Treatment and Procedure



Simulated Wood Cart

For a less industrial feel choose the Simulated Wood cart to match the décor of your facility.

UTWKA-33669-WOD
Lightweight Aluminum Unicart with wood veneer

Overall Dimensions: 29"W x 24.5"D x 42"H
 (Include casters, stabilizer base and ABS top)
 Drawer Dimensions: 22"W x 16.5"D
 Drawer Height: (2) 3"H, (2) 6"H, (1) 9"H
 Weight: 74 lbs.



Mobile Workstation


Access and capture clinical data and patient information at the point-of-care.



UTGSA-K33369-DKB
Lightweight Aluminum Unicart with blue gray drawers and key lock

Overall Dimensions: 29"W x 24.5"D x 42"H
 (Include casters, stabilizer base and ABS top)
 Drawer Dimensions: 22"W x 16.5"D
 Drawer Height: (1) 4" keyboard drawer
 (3) 3"H, (1) 6"H, (1) 9"H

The Mobile Workstation is shown with the following optional accessories:
 (1) KB-DRW – Keyboard drawer



Larger View – See Page 62

Phlebotomy/Specimen Collection Cart

FEATURES and BENEFITS

- These multipurpose carts are made out of a combination of steel, aluminum, and plastic to create a lightweight solution that is still very sturdy.
- Removable ABS plastic molded top.
- Includes gatelock for security.
- The 3 and 4 inch drawers act as tray holders for removable dividers. We include these customizable divider trays in each cart.
- Ergonomic push handles.
- 12 clear plastic tilt bins built into the front of the cart that conveniently show their contents for fast and easy access.
- Toss away waste into the pass-through sharps container that locks inside the cart or into the built-in tilt-out waste bin (see below).
- Also includes a glove box holder and side mounted paper holder.
- 5" Premium easy roll casters. (1) full swivel, (2) swivel w/ brake, (1) swivel w/ tracking.
- Stabilizer Base with high impact, non-marking plastic bumper to add stability and help keep cart from tipping.



MTWA-34696-WHT

Dimensions:

Overall: 27"W x 24"D x 45"H

Drawers:

- (1) 10.60" W x 16.5"D x 3"H
- (2) 10.60" W x 16.5"D x 4"H
- (3) 10.60" W x 16.5"D x 6"H
- (4) 10.60" W x 16.5"D x 9"H
- (5) 22"W x 16.5"D x 6"H

MSWA-3469-WHT

Dimensions:

Overall: 27"W x 24"D x 37.25"H

Drawers:

- (1) 10.60" W x 16.5"D x 3"H
- (2) 10.60" W x 16.5"D x 4"H
- (3) 10.60" W x 16.5"D x 6"H
- (4) 10.60" W x 16.5"D x 9"H




Unique lockable area to hold sharps container.



Built in drawer dividers for easy removal and cleaning.



Tilt out waste bin.



Paper holder & glove box holder.

Economy Carts by Waterloo




Features & Benefits

- Antimicrobial properties built-in to protect portions of the product.
- Centralized key locking standard. Breakaway locking kit available (Part # A15902A1)
- Drawers with ball bearing slides designed to push in and lock even after the cart is secured.
- New Ergonomic Push Handles for ease while pushing.
- Removable, easy to clean ABS plastic molded top.
- 5" Premium easy roll swivel casters (2 brake, 1 tracking, 1 full swivel).
- Stabilizer base with non-marking high impact plastic bumper to add stability.
- Integrated drawer handles run the length of the drawer for easy access.
- Additional accessories can be mounted to these carts, although this may require drilling holes.

See pages 32-39



ESGKU-33366-LTG



Enhance your look by adding any of our trim colors!

- TRIM-ECP - Pink
- TRIM-ECG - Green
- TRIM-ECB - Blue
- TRIM-ECV - Violet

Now available in compact size
20.5"W x 18"D x 36"H



Key Lock

Economy Carts



ESGKU-399-LTG

ESGKU-399-LTG Three-Drawer, Steel Economy Unicart

Dimensions
(include casters, stabilizer base and top)

Overall: 23"W x 24.5"D x 36"H

Drawers:

(1).....16.25"W x 16.5"D x 3"H

(2-3).....16.25"W x 16.5"D x 9"H

Weight: 78 lbs

ESGKU-3369-LTG Four-Drawer, Steel Economy Unicart

Dimensions
(include casters, stabilizer base and top)

Overall: 23"W x 24.5"D x 36"H

Drawers:

(1-2).....16.25"W x 16.5"D x 3"H

(3).....16.25"W x 16.5"D x 6"H


(4).....16.25"W x 16.5"D x 9"H

Weight: 84 lbs



ESGKU-3369-LTG

Economy Carts



WMT-4

We offer several styles of drawer dividers to help you organize your space.

WMT-4 Full drawer tray with (9) ampule dividers, ampule tray, rail and (4) adjustable dividers. Fits (1) per drawer. 16"W x 16."D x 2.3"H



WMT-5


WMT-5 Full drawer tray with (2) rails and (4) adjustable dividers. Fits (1) per drawer. 16"W x 16."D x 2.3"H

Elite Aluminum Uni-Cart Configurations

Waterloo Healthcare recognizes that there are unique and specific needs in most healthcare settings. Our catalog illustrates the most popular cart configurations, however we do not want to limit your options. Design your own configuration using the parameters outlined below.

Cart Size Options

1. Choose a cart size (Short, Medium, Tall, X-Tall)
2. Select any combination of 3", 6", 9" or 12" drawers that fit within your cart's configurable drawer space.
3. Select drawer front colors based on color options below.



Dimensions include new bumper, stabilizer base and casters.

Build Your Part Number, It's as Easy as 1,2,3...

EACH ALPHA OR NUMERIC DIGIT HAS A MEANING: example (UTGKA-33669-RED)

Section 1 letters

1st "U" = unicast
 2nd "S" = short, "M" = medium, "T" = tall
 or "X" = X-tall exterior
 3rd "B" = blue, "G" = lt. gray, "R" = red exterior
 "W" = wood, for Tall exteriors only
 4th "K" = key, "L" = lever, "S" = simplex,
 "E" = electronic, "P" = proximity
 5th "A" = aluminum

All carts shipped assembled.

Section 2 numbers

1st Height of 1st drawer
 2nd Height of 2nd drawer
 3rd Height of 3rd drawer
 4th Height of 4th drawer, etc.
 Totals for Section 2 must equal 21 inches for short cart;
 27 inches for tall cart; 30 inches for X-tall cart.


Section 3 colors (3 letters signify color choice for drawer fronts. Colors to choose from are listed below.)

| | | |
|------------------------|----------------------|-----------------|
| Charcoal.....GRA | NaturalNTL | VioletVIL |
| Dark BlueDKB | OrangeORG | WhiteWHT |
| Electric Blue.....ELB | Pink.....PNK | Wine.....WNE |
| Fairway GreenFWG | RedRED | YellowYEL |
| Light BlueLTB | Spring Green.....SPG | Wood.....WOD |
| Light Gray.....LTG | TealTLG | |


Elite Aluminum Med Junior Configuration

Cart Size Options


1. Choose a cart size (Short, Tall, X-Tall)
2. Select any combination of 3", 6", 9" or 12" drawers that fit within your cart's configurable drawer space.
3. Select drawer front colors based on color options below.




Med Jr. X-Tall
23"W x 24.5"D x 45"H



Med Jr. Tall
23"W x 24.5"D x 42"H



Med Jr. Medium
23"W x 24.5"D x 39"H



Med Jr. Short
23"W x 24.5"D x 36"H

30" Configurable Drawer Space 27" Configurable Drawer Space 24" Configurable Drawer Space 21" Configurable Drawer Space

Overall dimensions include new bumper and stabilizer base.

Build Your Part Number, It's as Easy as 1,2,3...

EACH ALPHA OR NUMERIC DIGIT HAS A MEANING: example (JSWKU-3369-FWG)

Section 1 letters

- 1st "J" = junior
- 2nd "S" = short, "M" = medium, "T" = tall
or "X" = X-tall exterior
- 3rd "G" = lt. gray, "W" = white exterior
- 4th "K" = key, "L" = lever, "S" = simplex,
"E" = electronic, "P" = proximity
- 5th "A" = aluminum

All carts are shipped assembled.

Section 2 numbers

- 1st Height of 1st drawer
- 2nd Height of 2nd drawer
- 3rd Height of 3rd drawer
- 4th Height of 4th drawer, etc.


Totals for Section 2 must equal 21 inches for short cart;
27 inches for tall cart; 30 inches for x-tall cart.

Section 3 colors (3 letters signify color choice for drawer fronts. Colors to choose from are listed below.)

| | | |
|------------------------|----------------------|-----------------|
| Charcoal.....GRA | NaturalNTL | VioletVIL |
| Dark BlueDKB | OrangeORG | WhiteWHT |
| Electric Blue.....ELB | Pink.....PNK | Wine.....WNE |
| Fairway GreenFWG | RedRED | YellowYEL |
| Light BlueLTB | Spring Green.....SPG | |
| Light Gray.....LTG | TealTLG | |


Traditional Steel Uni-Cart Configurations

Just like our NEW Elite Lightweight Aluminum carts, the Traditional Steel Uni-Carts are easily customizable to fit any healthcare facilities needs. With the following charts and instructions you can design a cart that will be unique and function specific. With so many options and no customization fee why go anywhere else?




X-Tall
29"W x 24.5"D x 46.5"H

27" Configurable Drawer Space
4" top drawer required.



Tall
29"W x 24.5"D x 43.5"H

24" Configurable Drawer Space
4" top drawer required.



Short
29"W x 24.5"D x 37.5"H

18" Configurable Drawer Space
4" top drawer required.

The top 4" drawer is required for our internal, tamper-resistant locking mechanism. Multiple 4" drawers are no longer available in any configuration. Dimensions include new bumper, stabilizer base and 5" casters.

Build Your Part Number, It's as Easy as 1,2,3...

EACH ALPHA OR NUMERIC DIGIT HAS A MEANING: example (UTGKU-433369-PNK)

Section 1 letters

- 1st "U" = uncart
- 2nd "S" = short or "T" = tall or "X" = X-tall exterior
- 3rd "B" = blue, "G" = lt. gray, "R" = red exterior
- 4th "K" = key, "L" = lever, "S" = simplex,
- 5th "U" = steel

All carts are shipped assembled.

Section 2 numbers

- 1st height of 1st drawer (A 4" drawer is standard here)
- 2nd Height of 2nd drawer
- 3rd Height of 3rd drawer
- 4th Height of 4th drawer, etc.
- Last digit of this section can be a "P" for a 9" panel.

Totals for Section 2 must equal 22 inches for short cart; 28 inches for tall cart; 31 inches for X-tall cart.


Section 3 colors (3 letters signify color choice for drawer fronts. Colors to choose from are listed below.)

| | | |
|------------------------|----------------------|-----------------|
| Charcoal.....GRA | NaturalNTL | VioletVIL |
| Dark BlueDKB | OrangeORG | WhiteWHT |
| Electric Blue.....ELB | PinkPNK | Wine.....WNE |
| Fairway GreenFWG | RedRED | YellowYEL |
| Light BlueLTB | Spring Green.....SPG | |
| Light Gray.....LTG | TealTLG | |

Traditional Steel Med Junior Configurations


Cart Size Options

1. Choose a cart size (Short, Tall, X-Tall)
2. Select any combination of 3", 6", 9" or 12" drawers that fit within your cart's configurable drawer space.
3. Select drawer front colors based on color options below.




Med Jr. X-Tall
23"W x 24.5"D x 45"H

30" Configurable Drawer Space



Med Jr. Tall
23"W x 24.5"D x 42"H

27" Configurable Drawer Space



Med Jr. Short
23"W x 24.5"D x 36"H

21" Configurable Drawer Space

Overall dimensions include new bumper and stabilizer base.

Build Your Part Number, It's as Easy as 1,2,3...

EACH ALPHA OR NUMERIC DIGIT HAS A MEANING: example (JSWKU-3369-FWG)

Section 1 letters

1st "J" = junior
 2nd "S" = short or "T" = tall or "X" = X-tall exterior
 3rd "G" = lt. gray, "W" = white exterior
 4th "K" = key, "L" = lever, "S" = simplex,
 "E" = electronic, "P" = proximity
 5th "U" = steel

All carts are shipped assembled.

Section 2 numbers

1st Height of 1st drawer
 2nd Height of 2nd drawer
 3rd Height of 3rd drawer
 4th Height of 4th drawer, etc.

*Totals for Section 2 must equal 21 inches for short cart;
 27 inches for tall cart; 30 inches for x-tall cart.*

Section 3 colors (3 letters signify color choice for drawer fronts. Colors to choose from are listed below.)

| | | |
|-------------------------|------------------------|------------------|
| Charcoal GRA | Natural NTL | Violet VIL |
| Dark Blue DKB | Orange ORG | White WHT |
| Electric Blue..... ELB | Pink PNK | Wine..... WNE |
| Fairway Green FWG | Red RED | Yellow YEL |
| Light Blue LTB | Spring Green SPG | |
| Light Gray LTG | Teal TLG | |



CBB-6



CBB-8



CBB-ER



OTB-4



OTB-G

UNI-CART® Accessories

CBB-6 Lightweight plastic cardiac board and metal mounting brackets for back side of Steel Unicart. Please use CBB-7 for Aluminum Unicarts, CBB-PC345 for PC343 & PC345. 24"W x 1"D x 18"H

CBB-8 Lightweight plastic cardiac board and side mounting brackets for back side of Aluminum Unicarts or Med Jr.

CBB-ER Lightweight plastic cardiac board mounting bracket and board for back side of the ER2000 cart.

OTB-4 Plastic oxygen tank brackets with adjustable velcro straps. For all Unicarts or Med Jrs. For PC343 or 345 please use part number OTB-PC345

OTB-G Green metal oxygen tank holder. Can be mounted to all Unicarts or Med Jrs. 4.75"W x 30"H

O-2 Oxygen tank bracket with adjustable velcro strap for ER2000 carts. Included with purchase of ER2000.

SB-1 Sharps container mounting bracket with adjustable velcro strap. Mounts to side of all Unicarts, Med Jr, ER-2000, Medtrix, DMC.

SB-2L Locking sharps container mounting bracket. Mounts to side of all Unicarts, Med Jr, Medtrix, DMC. Designed for 5 quart container. Also available in non-locking version under SB-2

UH-2 Two utility hooks. Mounts to existing holes on either side or back of cart.

ABH-2 Set of two self-adhesive chrome hooks that mount easily on any cart.



O-2



SB-1



SB-2L



UH-2



ABH-2

UNI-CART® Accessories



IV-2

IV-2 Chrome IV pole and mounting bracket. Mounts on Unicart, Med Jr, Medtrix, DMC. For ER2000 please use part number IV-3.



CR-1

CR-1 Coat rack with hooks. Made from lightweight aluminum. Mounts on Unicart, Med Jr, Medtrix, or DMC carts.



BH-3

BH-3 The NEW glove box holders are available in 3 sizes to fit any facilities needs. Order BH-1 for one box (10"H x 4"D x 6"W), BH-2 for two box size (10"H x 4"D x 12"W), BH-3 for three box size (10"H x 4"D x 18"W)

GMH-1 Glove & Mask box holder. Made from lightweight aluminum. Adjusts to many sizes of boxes. Mounts to railing on cart.

CMB-3 Mounting brackets for mounting cassette products over Steel Unicart. For cassettes, please see CASS-4S (4-tier), CASS-3S (3-tier), CASS-2S (2-tier).



GMH-1

CAST-1 Side mounting bucket holder and towel bar. Holds (2) 165oz cast buckets. 17.5"W x 6.5"D x 8.8"H

BSK-1 Small (6 quart) plastic waste basket. Requires accessory rail. 9"W x 6.25"D x 10"H

BSK-2 Large (20 quart) plastic waste basket. Requires accessory rail. 13"W x 8.5"D x 15"H

KB-DRW Keyboard drawer for Steel Unicarts may be purchased as an option in place of a 4"H drawer. Aluminum carts purchase KA-DRW as an option in place of a 3"H drawer.



CAST-1



BSK-1



BSK-2



KB-DRW



TB-3, TB-4, TB-5, TB-6




UR-1, UR-2, UR-3



CHT-1



WNH-1



WSB-1

UNI-CART® Accessories

TB-3, TB-4, TB-5, TB-6 Clear plastic tilt bins with white exterior. Available in 3, 4, 5, or 6 bin units. Mounts to base rail system (BRS-2, BRS-3 OR BRS-JR)

UR-1, UR-2, UR-3 Utility receptacles, available in three sizes, UR-1 (0.85 liters), UR-2 (1.5 liters), UR-3 (2.8 liters). Requires accessory rail.

CHT-1 Gray steel chart holder. Mounts to accessory rail. 14.25"W x 3.75"D x 10"H

WNH-1 Wire notebook holder. Mounts to accessory railing. 12.5"W x 2.5"D x 10"H

WSB-1 Extra-large wire basket. Mounts to accessory railing. 15.5"W x 9.5"D x 5.5"H

WB-1L Large wire basket. Mounts to accessory railing. 15.5"W x 5.5"D x 3.5"H

WB-1S Small wire basket. Mounts to accessory railing. 8.5"W x 5.5"D x 3.5"H

ORG-1L Large clear plastic organizer. Includes clear slide-in dividers. Mounts to accessory railing. 18"W x 5"D x 3"H

ORG-1S Small clear plastic organizer. Includes clear slide-in dividers. Mounts to accessory railing. 9"W x 5"D x 3"H

TD-1 Light gray tape dispenser with suction cups to attach to the top of a cart.



WB-1L




WB-1S



ORG-1L



ORG-1S




TD-1

UNI-CART® Accessories



BRS-2



R-3



SHLF-2



SHLF-110



ER-SHLF

BRS-2 Horizontal rail system for hanging accessories and mounting shelves on steel Unicarts. Includes two vertical rails and three horizontal rails. Med Jr. carts purchase BRS-JR. Aluminum carts purchase BRS-3.

R-3 Wrap-around aluminum railing for Steel Unicarts. May be used to hang accessories. Med Jr. carts purchase R-JR. Aluminum carts purchase R-4.

SHLF-2 Flat shelf that mounts to the horizontal rail system (BRS-2). Aluminum carts purchase SHLF-3 and BRS-3. 24"W x 10"D x 3/4"H

SHLF-110 Two-tier shelving unit that mounts above Steel Unicarts. Order **SHLF-110A** for aluminum carts. 26.5"W x 15.5"D x 1"H

ER-SHLF Full swivel monitor shelf capable of locking in place for Steel carts. Includes an adjustable strap to secure items in place. Maximum weight capacity 25 lbs. Instant Access series purchase ER-SHLF2000. Aluminum carts purchase ER-SHLFA. 13"W x 17"D x 12"H

SH-4 Our durable steel scope holders secure easily to the wall or our Unicarts. Available in two sizes, two scope holder (SH-2) and four scope holder (SH-4). Locking Plexiglas doors also included for added security.

SSG-2C Collapsible side shelf for Steel Unicart and Med Jr. Hangs on either side of cart. Surface area is 17.5"W x 15"D

SS-M Suction shelf for Unicart or Med Jr. 18"W x 8"D x 3"

CH-3 Catheter holder for Unicart or Med Jr. 7.5"W x 2"D x 23"H


SUT-1 This shelving unit includes (2) vertical rails, (1) top horizontal rail and (2) suture shelves. Aluminum carts require the purchase of AMP-2 for installation of this item. Each shelf measures 21.5"W x 5.5"D x 3.25"H




SH-4



SSG-2C



SS-M




CH-3




SUT-1



AMP-1



AMP-2



SR-1



CB-1



OS-1

UNI-CART® Accessories

AMP-1 Accessory mounting panel for Lightweight Aluminum Unicarts. Mounts to **SIDE** to hold additional accessories.

AMP-2 Accessory mounting panel for Lightweight Aluminum Unicarts. Mounts to **BACK** to hold additional accessories.

SR-1 Accessory rail for side attachment of accessories to the Steel Unicarts or Med Jr. Aluminum carts purchase AMP-1.

CB-1 Accessory rail for back attachment of accessories to the Steel Unicarts only. Aluminum carts purchase AMP-2.

OS-1 Hospital grade outlet strip. 15 foot long cord with 6 outlets.

CW-1 Cord wrap. Usually used with outlet strip (OS-1) sold separately.

FT-1 Fluid/equipment storage tray system for Steel, Med Jr or Medtrx/Halftx carts. Includes three trays. Aluminum carts purchase FTA-1. Each tray measures 17.52"W x 4.04"D x 4.5"H

T-2 ABS plastic top for Steel Unicarts. Easily slips on and off for cleaning. 27" x 18.5". Med Jr carts purchase T-JR. Aluminum carts Purchase T-3.

CP-1 Replacement caster package. 5" premium easy roll swivel casters. Includes: (2) Brake w/swivel, (1) tracking w/swivel and (1) full swivel caster.

BAG-COVER5 Plastic cart cover for Unicarts. Dimensions are 33" x 23" x 42".


LH-1 Label holders, self-adhesive. 1"H x 3"W, package of 100. Holders are loaded from the side.



FT-1



T-2



CP-1



BAG-COVER5



LH-1

UNI-CART® Security Options



SD-3

SD-3 Locking security cabinet for use of transportation of narcotics or permanently mounted inside drawer. Includes two keys. 15.6"W x 11.4"D x 2.5"H



PROX-1

ELECTRONIC & PROXIMITY LOCKING SYSTEMS

Increase the security of your carts from the start! Electronic locking systems available on our Med Junior and Aluminum Unicarts. Upgrade with our proximity reader and use your badge for access! Call 1-800-833-4419 to learn more.



PROX-2

LOCKMANAGER Management software for Electronic and Proximity locking systems.

NARC-BOX Narcotic Box with security lid to secure medication in any full size Steel or Aluminum Unicart drawer 9"H. Includes lock with two keys.



PS-100

PROX-1 Proximity card. 2"W x 3"H

PROX-2 Small, adhesive proximity tag. 1" diameter.

PS-100 Plastic breakaway seals. Can be used with lever locking systems, breakaway gate locks and EMER-BOX divider trays.



ELECTRONIC & PROXIMITY LOCKING SYSTEMS



NARC-BOX

Uni-Cart Locking Options

Gate Locking Options



GLB

ELITE GATE LOCKS are for use on any of our lightweight aluminum carts. Can be used for securing narcotics and liquids or as a visible inventory management system.

STEEL GATE LOCKS are for use on our traditional steel carts. Assembly required. If placing an order for a cart just let us know and we will attach it to the cart prior to shipping!

Call 1-800-833-4419 and we will help you decide what lock configuration will work best for your needs.




GLK

Dividers/Trays



CP-DIV-TRAY



DIV-CC2



WDC-2025



WDC-2050



EMER-BOX

CP-DIV-TRAY Clear plastic divider trays (two per kit). Includes dividers. Fits in all drawers except half-sized Unicart or Medtrix drawers. 8.5" x 14" x 2"

DIV-CC2 & DIV-CC4

Cut and clip divider systems. Simply score the plastic strips and snap them to create your own divider system. Fits in all drawers. Set includes (6) 22" long white plastic strips and 8 clips. DIV-CC2 include 2" H strips and DIV-CC4 include 4"H strips. Additional clips available in sets of 6 (DIV-CLIP, shown at left in DIV-CC2 image).

WDC-2025 Configurable gray plastic divider tray. 2.5 inches deep, includes dividers. Fits two per drawer, and typically used in 3" and 4" deep drawers. 9.25"W x 15"D x 2.5"H

WDC-2050 Configurable gray plastic divider tray. 4.75 inches deep, includes dividers. Fits two per drawer, and typically used in 6", 9", and 12" deep drawers. 9.25"W x 15"D x 4.75"H

EMER-BOX White plastic divider box with clear plastic lid. Can be sealed with plastic seal. 21.75"W x 15.75"D x 2.5"H.

TT-1A White plastic tray with configured bins. Can be reconfigured easily. Fits in all full sized Unicart drawers. 22"W x 15"D x 2.5"H

TT-1B White plastic tray with configured bins. Can be reconfigured easily. Fits in all full sized Unicart drawers. 22"W x 15"D x 2.5"H

DIV-DRW12 Adjustable divider system with a row of ampule holders. 21.5"W x 15.25"D. To purchase this divider system for the ER2000, instant access cart, please order DIV-DRW12ER.

DIV-DRW34 Adjustable divider system with large divided sections. 21.5"W x 15.25"D. To purchase this divider system for the ER2000, instant access cart, please order DIV-DRW34ER.


DIV-DRW5 Adjustable divider system with large divided sections. Fits in 6, 9, 12 inch deep full sized Unicart drawers. 21.5"W x 15.25"D x 5"H. To purchase this divider system for the ER2000, instant access cart please order DIV-DRW5ER.



TT-1A



TT-1B



DIV-DRW12



DIV-DRW34



DIV-DRW5

Dividers/Trays



DIVTRAY-AMP

DIVTRAY-AMP Ampule holder divider system with white plastic tray. Fits in all full sized Unicart drawers. 21.5"W x 15.25"D x 2.3"H



DIVTRAY-ANES

DIVTRAY-ANES Anesthesia divider system with white plastic tray. Fits in all full sized Unicart drawers. 21.5"W x 15.25"D x 2.3"H




DIVTRAY-STD

DIVTRAY-STD Standard divider system with white plastic tray. Fits in all full sized Unicart drawers. 21.5"W x 15.25"D x 2.3"H. Also available in 5"H, purchase DIVTRAY-STD5.



WMT-3

WMT-3 Complete ampule holder divider system with white plastic tray. Includes (3) rows of 17 ampule dividers. 21.5"W x 15.25"D x 2.3"H




WMT-4

Slim Med Junior Trays

WMT-4 Full drawer tray with (12) ampule dividers, ampule tray, rail and (4) adjustable dividers. Fits (1) per drawer. 16"W x 16"D x 2.3"H

WMT-5 Full drawer tray with (2) rails and (4) adjustable dividers. Fits (1) per drawer. 16"W x 16"D x 2.3"H



WMT-5




DIV-UC6

DIV-UC6 & DIV-UC9 Front to back divider system for 6" or 9" Unicart drawers (Not for Med Jr. drawers). Includes 3 front to back panels to create up to 4 compartments in one drawer.



ER-TRAY3 & ER-TRAY6

ER-TRAY3 & ER-TRAY6 Clear trays for Unicart & Instant access series (ER2000). For dividing systems purchase the DIV-DRW12ER, DIV-DRW34ER or DIV-DRW5ER



BAG-3 & BAG-6

BAG-3 & BAG-6 Package of 20 sealable perforated bags that fit over divider trays. Dimensions: BAG-3: 23"W x 17"D x 5.5"H BAG-6: 23"W x 17"D x 7"H

DIV-10L Long clear plastic additional dividers for DIV-DRW12, DIV-DRW34, DIVTRAY-STD, DIVTRAY-AMP, DIVTRAY-ANES.


DIV-10S Short clear plastic additional dividers for DIV-DRW12, DIV-DRW34, DIVTRAY-STD, DIVTRAY-AMP, DIVTRAY-ANES.



DIV-10L & DIV-10S

Continued on next page

Dividers/Trays



DIV-AU3




DIV-AU6

DIV-AU3 Divider kit for 3" Aluminum Unicart Drawers. Includes one front-to-back and one side-to-side divider.

DIV-AU6 Divider kit for 6", 9", 12" Aluminum Unicart Drawers. Includes one front-to-back and one side-to-side divider.

DIV-AJ3 Divider kit for 3" Aluminum Med Jr Drawers. Includes one front-to-back and one side-to-side divider.

DIV-AJ6 Divider kit for 6", 9", 12" Aluminum Med Jr Drawers. Includes one front-to-back and one side-to-side divider.



DIV-AJ3



DIV-AJ6

Accessory Packages

AIR-PKG Difficult Airway (breathing) accessory package for Traditional Steel Unicarts includes: Side shelf (SSG-2C) Vertical rail system (BRS-2), Shelf for rail system (SHLF-2) Small organizer, clear with lid (ORG-1S), Two oxygen tank brackets (OTB-4).

AIR-PKGA Difficult Airway (breathing) accessory package for Elite Aluminum Unicarts includes: Vertical Rail System (BRS-3), Shelf for rail system (SHLF-3) Small organizer, clear with lid (ORG-1S), Two oxygen tank brackets (OTB-4).

ANES-PKG Standard Anesthesia accessory package for Traditional Steel Unicarts. See Page 17 (cart sold separately)

ANES-PKGA Standard Anesthesia accessory package for Elite Aluminum Unicarts. See Page 16 (cart sold separately)

ANES-PKG-DLX Deluxe Anesthesia accessory package for Traditional Steel Unicarts. See Page 17 (cart sold separately)

ANES-PKGA-DLX Deluxe Anesthesia package for Elite Aluminum Unicarts. See Page 16 (cart sold separately)

EMER-PKGM Emergency / Crash accessory package for Traditional Steel Unicarts. See Page 12 (cart sold separately)

EMER-PKGA Emergency / Crash accessory package for Elite Aluminum Unicarts. See Page 12 (cart sold separately)

EMER-PKGM-DLX Deluxe Emergency / Crash accessory package for Traditional Steel Unicarts. See Page 11 (cart sold separately)

EMER-PKGA-DLX Deluxe Emergency / Crash accessory package for Elite Aluminum Unicarts. See Page 11 (cart sold separately)

INF-PKGM Infection / Isolation accessory package for Traditional Steel Unicarts includes: Coat rack (CR-1), 2 each Utility hooks (UH-2) and Railing and supports (R-3).

INF-PKGA Infection / Isolation accessory package for Elite Aluminum Unicarts includes: Coat rack (CR-1), 2 each Utility hooks (UH-2) and Railing and supports (R-4).

IV-PKG IV accessory package for Traditional Steel Unicarts includes: IV pole and bracket (IV-2), Sharps bracket (SB-1), Tape dispenser (TD-1), Base rail system (BRS-2) and 4 compartment tilt bins (TB-4).

PEDS-PKG Pediatric accessory package for Traditional Steel Unicarts includes: Base rail system (BRS-2), Waste receptacle (BSK-1), IV pole and bracket (IV-2), Top shelf (SHLF-2), Side rail (2 each SR-1), Chart holder (CHT-1) Collapsible side shelf (SSG-2C), 5 compartment tilt bins (TB-5), Railing and supports (R-3).

PEDS-PKGA Pediatric accessory package for Elite Aluminum Unicarts includes: Base rail system (BRS-3), Top Shelf (SHLF-3), Tilt bins (TB-5), Waste container (BSK-2), Railing and supports (R-4), Chart holder (CHT-1), IV pole and bracket (IV-2), Plastic breakaway seals (PS-100), (9) 3" gatelock bar (GLB-3).

SUT-PKG Suture accessory package for Traditional Steel Unicarts includes: Wrap-around railing (R-3), (2) Suture shelves (SUT-1), and two Utility hooks (UH-2).

Economy and Portable Storage

PC-343



Waterloo also offers carts for facilities that require an economical solution due to budget demands. Our Economy Carts are ideal for sports medicine, veterinary supply carts, or any setting requiring secure storage. Included breakaway bar can be used down the front of cart with our PS-100 seals or a personal padlock.

We offer accessories that can be mounted with minimal installation.

- CBB-PC345** Cardiac Board and Brackets
- OTB-PC345** Oxygen Tank Brackets
- ABH-2** Adhesive Utility Hooks

Other accessories found on pages 32-37 can be mounted on these carts, although this may require drilling holes.

PC-345



Dimensions (includes casters, bumper and top)

- Overall34"W x 20"D x 34.8"H
- Drawers
- (1)22"W x 16.8"D x 2"H
- (2-3).....22"W x 16.8"D x 4"H
- Panel.....22"W x 16"D x 11.5"H

Dimensions (includes casters, bumper and top)

- Overall34"W x 20"D x 34.8"H
- Drawers
- (1)22"W x 16.8"D x 2"H
- (2-3).....22"W x 16.8"D x 4"H
- (4-5).....22"W x 16.8"D x 6"H

You can always be prepared with our portable storage by your side. These lightweight totes can hold essential supplies for immediate use. Available with caddy to assist in transporting.



EMB-1C Stackable Three Compartment Tote

Caddy included for easy transport w/tub, (2) trays and lid. Accepts seals or padlock for security. 13.5"W x 21"D x 17"H

MINI-BOX Compact Three Compartment Tote

Allows you to have supplies on hand at any site. Accepts seals or padlocks for security. 16.25"W x 9"D x 12.5"H

Narcotics Cabinets

When you have a need to keep security tight on medications and supplies, Waterloo Healthcare offers the perfect solution with its full line of **narcotic cabinets** with single or dual security locks. These cabinets are designed for use in office areas, clinics, or to store supplies for specific procedures in any area where patients are receiving care.



NCS1

Please note. Cabinets do NOT come keyed alike. If you prefer the cabinets keyed alike, make note on your purchase order.

Small Cabinets

Overall Dimensions: 12"W x 8"D x 15"H

NCS1 Single Door/Single Lock
Includes three adjustable or removable shelves.

NCS2 Double Doors/Double Lock
Includes three adjustable or removable shelves, double locks keyed differently.



NCL-D1L1



NCL-D1L2



NCS2

Large Cabinets

Overall Dimensions: 24"W x 10"D x 30"H

NCL-D1L1C Single Door/Combination Lock.
Includes three shelves: one permanent and two adjustable or removable. Combination lock with key override.

NCL-D1L1 Single Door/Single Lock
Includes three shelves: one permanent and two adjustable or removable.

NCL-D1L2 Single Door/Double Lock
Includes three shelves: one permanent and two adjustable or removable, double locks keyed differently.

NCL-D2L2 Double Door/Double Lock
Includes two adjustable or removable shelves, double locks keyed differently. 18"W x 8"D x 30"H



NCL-D2L2



NCL-D1L1C

Narcotics Cabinets

Security Drop Box

- Secure envelopes, keys, confidential paperwork or specimen samples for lab pickup.
- Items are deposited through the front opening for select personnel retrieval. A mail box type baffle deters the removal of the contents.
- Secured by two locks, keyed differently.
- Overall Dimensions: 14.5"W x 20"H x 17"D



SDB-1

Medication Cart Systems

Waterloo Healthcare medication carts (trolleys) offer quality...security...and...flexibility.



MEDTRX 1.5

Dimensions (includes casters)
Overall exterior size 39"W x 20"D x 43"H



MEDTRX-3AC

*Shown with optional narcotics drawer lock (NL-MT)
Pictured as part #MTK14B771/TFTFT*

MEDTRX 1.0

Dimensions (includes casters)
Overall exterior size 27"W x 20"D x 43"H



HALFTRX-2LTC

*Pictured as part #HTK1-999
Holds 450 pill cards.*

CASSETTE UNICART



CASSETTE UNICART

*Pictured above as part
#UTGKU-C3336-LTG*

DMC-3 (includes side module)

Dimensions (Includes stabilizer base, side handle and casters)
Overall32"W x 24.5"D x 43.5"H



DMC3-4S10

PHOENIX LTC CARTS

Overall exterior size: 38"W x 24"D x 43"H



750PC
Holds 750 pill cards.

Medication Cart

The NEW Aluminum Wood Simulated Cart from Waterloo Healthcare



Features & Benefits

- Antimicrobial properties built-in to protect portions of the product
- Elite Aluminum offers easier maneuverability and safety with a lightweight aluminum construction being 30% lighter than standard steel carts
- ABS plastic molded top.
- New dual pull-out shelves provide additional work surface on both sides of the cart.
- Adjustable dual push handles.
- Centralized locking with 5 available options.
- Drawers with ball bearing slides designed to close and lock even after the cart has been secured.
- Drawers include built-in, interlocking dividers corresponding to the drawer height.
- 5" Premium easy roll swivel casters (2 brake, 1 tracking, 1 full swivel).
- Stabilizer base with high impact plastic bumper to add stability.
- Latex Free.
- Available in one size:
Tall: (27" drawer height)
29" W x 24.5" D x 42" H
and 12-inch deep drawers.



Simulated Wood Cart

For a less industrial feel choose the Simulated Wood cart to match the décor of your facility.

Customize to use for unit dose boxes or use as a pill card cart with the ability to configure with 3 inch, 6 inch, 9 inch or 12 inch deep drawers.

UTWKA-33669-WOD

Lightweight Aluminum Unicart with wood veneer.

Overall Dimensions: 29"W x 24.5"D x 42"H
(Include casters, stabilizer base and ABS top)
Drawer Dimensions: 22"W x 16.5"D
Drawer Height: (2) 3"H, (2) 6"H, (1) 9"H

Weight: 74 lbs.

Aluminum Locking Options
Page 6

Aluminum Accessory Options
Pages 32-37

Self Configuration - Tall Only
Page 29

UNIT DOSE CART

UXGSA-333333336-LTG

(Cart shown with optional pushbutton locking system).

Dimensions

(includes casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 45"H

Drawers: 22"W x 16.5"D

Drawer Height: (8) 3"H, (1) 6"H

Weight: 85 lbs

5 available locking systems. See Page 6.



UXGSA-333333336-LTG



UXGKA-3999-LTG

PILL CARD CART

UXGKA-3999-LTG

Dimensions

(includes casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 45"H

Drawers: 22"W x 16.5"D

Drawer Height: (1) 3"H, (3) 9"H

Weight: 85 lbs

5 available locking systems. See Page 6.

Holds 450 pill cards.



DIV-UC9

DIV-UC9 Divide the drawer into rows to provide organization for pill cards or fluids. For additional Aluminum cart accessories see page 32 - 39.



NARC-BOX

NARC-BOX Secure medications in the 9"H aluminum or steel Unicart drawers.

Larger View - See Page 62



GRA-Charcoal



DKB-Dark Blue




ELB-Electric Blue



FWG-Fairway Green



LTB-Light Blue



LTG-Light Gray



NTL-Natural



ORG-Orange



PNK-Pink



RED-Red




SPG-Spring Green



TLG-Teal



VIL-Violet



WHT-White



WNE-Wine



YEL-Yellow

Medication Cart Highlights

Waterloo Healthcare medication carts offer quality, security and flexibility. With two locking options and flexible custom configuration, the Medtrix-series carts will fit the needs of any facility. See page 57 for customer configuration options.

Our Medtrix series is available in two sizes. Our full size cart offers half size drawers on the left for storage and supplies and space for up to 50 patient bins or full size drawers on the right. We also have a half size cart (Halftrix) for use where space is limited. This cart can hold full size storage drawers, patient bins or a mix of both.

Locking options to choose from:



Tamper Resistant Best® Key Locking System



Simplex Pushbutton With Key Lock Override

PATIENT CASSETTE CART SYSTEMS

Styled push handles

Tamper resistant Best® Lock

Open/closed warning label on both sides

Optional narcotic lock (not available on 3"H drawers)



MEDTRX-1AC

Pictured as part #MTK1-4B411/6TTTT

PILL CARD CART SYSTEMS

Antimicrobial properties built-in to protect portions of the product.



Antimicrobial Surfaces



HALFTRX-1LTC

Pictured as part #HTK1-33311

Slide out shelf

ABS plastic removable molded top

Pill card drawers hold a total of 300 pill cards.

Wrap around bumper

5" premium easy roll casters (one swivel, one tracking, two w/locking toe brake)

*Quantity will vary depending on card thickness and use of dividers.

Medication Cart Highlights

When dual-sided accessibility is required, look no further, a time proven design from Waterloo Healthcare, our dual-sided medication carts have been popular worldwide for many years.

PATIENT CASSETTE CART SYSTEMS

Pull-out shelf DMC3-4S10



Both sides of cart have same drawer configuration.

- Cup holder
- Push handle
- Positive center catch – prevents drawers from accidentally passing through
- Lockable/portable cassette storage
- Multiple drawer sizes
- Stabilizer Base with Bumper Guard adds stability by lowering the center of gravity on the cart, keeping it from tipping
- 5” premium easy roll casters (one swivel, one tracking, two w/locking toe brake)

Antimicrobial properties built-in to protect portions of the product.




DMC-3 Dual-Sided Medication Cart
 Dimensions (includes casters)
 Overall32”W x 24.5”D x 43.5”H

SIDE VIEW



DRF-2

Add lockable security to storage area on side module. Ideal for liquid storage.



- Removable, replaceable ABS top
- Tamper resistant Best® Lock
- Storage shelves
- Dual access drawers
- Roller bearing slides
- Tilt waste receptacle

Medication Cart

MEDTRX & HALFTRX Patient Cassette Cart Systems

Select one of these configurations or build your own custom cart, see page 57. All carts shown with optional narcotics lock (NL-MT) and drawer trim (see page 56).



Antimicrobial Surfaces



MEDTRX-1AC

Pictured as part #MTK1-4B411/6TTTT



MEDTRX-2AC

Pictured as part #MTK1-4B411/TFTF6



MEDTRX-3AC

Pictured as part #MTK1-4B771/TFTFT



MEDTRX-4AC

Pictured as part #MTK1-4B771/33FTF6



MEDTRX-5AC

Pictured as part #MTK1-4B771/33TTTT



HALFTRX-1AC

pictured as part #HTK1-6BTTTT



HALFTRX-2AC

pictured as part #HTK1-33TTT6

Dual-Sided Medication Carts

Standard Features:

- Tamper resistant Best® locking mechanism maintains secured storage of all components.
- Bin sizes available: 3", 4", 6" wide bins available.
- Cassette shelves in a DMC-3 house two cassettes of the same size, back to back.
- Four 5" premium easy roll casters allow you to move from tile to carpeted surfaces with ease. Comes with two swivel w/ brake, one full swivel and one tracking caster
- Two shelves and built-in waste receptacle.
- Non-marking, high-impact stabilizer base.
- Slide out work surface with cover, ABS plastic top and cup dispenser are also standard with the DMC carts.
- Drawers can be opened from either side of cart. Ball-bearing drawer slides provide easy opening and are unaffected by weight of stored items. A positive center catch device prevents drawers from accidentally passing through when cart is being moved or when unlocked.

DMC3-4S10



Antimicrobial Surfaces



Side Module w/Shelf
Additional storage.



Waste Receptacle
Dispose of your waste easily with this built-in receptacle, conveniently located in side module.

Dual-Sided Med.

Phoenix LTC Pill Card & Unit Dose Box Carts



450PC
Holds 450 pill cards.



600PC
Holds 600 pill cards.



750PC
Holds 750 pill cards.



900PC
Holds 900 pill cards.



198UDB
Holds 198 unit dose boxes.



252UDB
Holds 252 unit dose boxes.



306UDB
Holds 306 unit dose boxes.



378UDB
Holds 378 unit dose boxes.

**For more product information or accessories,
please visit www.phoenixltc.com**

MEDTRX & HALFTRX Pill Card Cart Systems

Select one of these configurations or build your own custom cart, see page 57. Some carts shown with optional narcotics lock (NL-MT) and drawer trim (see page 56).



MEDTRX-1LTC
Pictured as part #MTK1-4B411/33311
 Holds 400 pill cards.



MEDTRX-2LTC
Pictured as part #MTK1-4B411/3611
 Holds 400 pill cards.



MEDTRX-3LTC
Pictured as part #MTK1-4B411/999
 Holds 550 pill cards.



MEDTRX-4LTC
Pictured as part #MTK1-999/33311
 Holds 450 pill cards.



MEDTRX-5LTC
Pictured as part #MTK1-4771/999
 Holds 500 pill cards.



MEDTRX-9LTC
Pictured as part #MTK1-999/999
 Holds 600 pill cards.



HALFTRX-1LTC
Pictured as part #HTK1-33311
 Holds 300 pill cards.




HALFTRX-2LTC
Pictured as part #HTK1-999
 Holds 450 pill cards.



Antimicrobial Surfaces


MEDTRX & HALFTRX




CASS-3MTC



CASS-5MTC



CASS-7MTC



BIN-3G



BIN-5G

Medtrix / Halftrix Cassettes and Bins

CASS-3MTC Two-tier cassette includes (10) 3" wide bins for MEDTRX series carts, exchange programs or individual use. Additional bins sold separately (BIN-3G) 16.5"W x 13.5"D x 6"H.

CASS-5MTC Two-tier cassette includes (6) 5" wide bins for MEDTRX series carts, exchange programs or individual use. Additional bins sold separately (BIN-5G) 16.5"W x 13.5"D x 6"H).

CASS-7MTC Two-tier cassette includes (4) 7" wide bins for MEDTRX series carts, exchange programs or individual use. Additional bins sold separately (BIN-7G) 16.5"W x 13.5"D x 6"H.

BIN-3G Patient bins for MEDTRX medication cassettes (10/ctn).
2 7/8"W x 10 3/4"D x 2 3/16"H

BIN-5G Patient bins for MEDTRX medication cassettes (6/ctn).
4 13/16"W x 10 3/4"D x 2 3/16"H


BIN-7G Patient bins for MEDTRX medication cassettes (4/ctn).
7 1/4"W x 10 3/4"D x 2 3/16"H

BD-3B Left to right bin dividers that fit 3" MEDTRX bins. Fits part number BIN-3G.


BD-5B Left to right bin dividers that fit 5" MEDTRX bins. Fits part number BIN-5G.

BD-7B Left to right bin dividers that fit 7" MEDTRX bins. Fits part number BIN-7G.


LBL-1B Bin labels for MEDTRX bins. 1000 labels. Size: 2.5"W x 0.75"H




BIN-7G



BD-3B



BD-5B



BD-7B



LBL-1B

DMC Cassettes and Bins



CASS-2S




CASS-3S



CASS-4S



MBG-3B



MBG-4B

CASS-2S Two-tier lockable cassette for dual-sided medication carts, exchange programs or individual use. Bins sold separately. Multiple orders can be "keyed alike" just use part number CASS-2S110

21"W x 10"D x 7.75"H

BINS SOLD SEPARATELY

CASS-3S Three-tier lockable cassette for dual-sided medication carts, exchange programs or individual use. Bins sold separately. Multiple orders can be "keyed alike" just use part number CASS-3S110

21"W x 10"D x 10.75"H

BINS SOLD SEPARATELY

CASS-4S Four-tier lockable cassette for dual-sided medication carts, exchange programs or individual use. Bins sold separately. Multiple orders can be "keyed alike" just use part number CASS-4S110

21"W x 10"D x 13.75"H

BINS SOLD SEPARATELY

MBG-3B 3" wide patient bins for DMC medication cart (24/ctn). 3"W x 9"D x 2.75"H. Dividers sold separately.

MBG-4B 4" wide patient bins for DMC medication cart (16/ctn). 4.5"W x 9"D x 2.75"H. Dividers sold separately.


MBG-6B 6" wide patient bins for DMC medication cart (12/ctn). 6"W x 9"D x 2.75"H. Dividers sold separately.

BDM-2 Left to right bin dividers that fit 4.5" DMC bins. Fits part number MBG-4B. Bins sold separately.

BDS-1 Left to right bin dividers that fit 3" DMC bins. Fits part number MBG-3B. Bins sold separately.

BDW-6 Left to right bin dividers that fit 6" DMC bins. Fits part number MBG-6B. Bins sold separately.

PBL-1 Bin labels for plastic bins in DMC carts (MBG-3B, 4B, 6B). 1 roll of 1000 labels for individual patient bins.



MBG-6B



BDM-2



BDS-1




BDW-6



PBL-1


DMC Cassettes



DIV-HTS



DIV-HTL



DIV-MTS



DIV-MTL

Medtrx / Halftrx Pill Card Dividers

DIV-HTS Divides half-size 4" - 7" drawers from front to back. Includes dividers and brackets (3.5"H).


DIV-HTL Divides half-size 9" or larger drawers from front to back. Includes dividers and brackets (8.5"H).

DIV-MTS Divides full-size 4" - 7" drawers from front to back. Includes dividers and brackets (3.5"H).

DIV-MTL Divides full-size 9" or larger drawers from front to back. Includes dividers and brackets (8.5"H).



LC-1C



DRF-2

Medtrx / Halftrx Medication Security

LC-1C Full or half width locking covers are available for added security and can be "keyed alike" or different from cart. Double hinge in the center allows access from either side. (Half size cover available for 2.75" and 4" drawers only.)

DRF-2 Add lockable security to storage area on side module. Ideal for liquid storage

NL-1.0C Half-sized lockable security drawer cover for MEDTRX series carts. Can be installed on 4" deep or larger.

NL-1.5C Full-sized lockable security drawer cover for MEDTRX series carts. Can be installed on 4" deep or larger.

NARC-BOX Narcotic Box with security lid to secure medication in any full size Steel or Aluminum Unicart drawer 9"H. Includes lock with two keys



NL-1.0C



NL-1.5C



NARC-BOX

Medtrix / Halftrix Accessories



PC1000 & PC3000

PC1000 1000 Pill crushing replacement pouches. Tablets are crushed in a durable latex-free pouch.

WL-3 Screws into wall and allows you to secure your medication cart to the wall.

CDA-2 Adjustable pop-up style cup dispenser. Holds 2 to 5 oz. cups. For use with all of our carts.

EWR-1 Wire trash container frame with lid and 10 bags. Replacement bags available (EWR-BAG) 12.5"W x 7"D x 19.5"H

BSK-1 Small (6 quart) plastic waste basket. Requires mounting rail (SR-1). 13"W x 6.3"D x 10"H

BSK-2 Large plastic waste basket. Requires mounting rail (SR-1). 13"W x 8.5"D x 15"H

SR-1 A side rail for hanging accessories on sides of carts.


SB-2L Sharps bracket includes key lock for added safety & security. Fits 5 quart sharps container.

SB-1 Sharps container mounting bracket with adjustable velcro strap. Mounts to side of all Unicarts, Med Jr, ER-2000, Medtrix, DMC.

LH-1 Label holders, self-adhesive. 1"H x 3"W, package of 100. Holders are loaded from the side.



BSK-2



SR-1



SB-2L



SB-1



LH-1



WL-3



CDA-2




EWR-1



BSK-1



OS-1



CW-1



UR-1, 2, & 3



IV-2



FT-1

Medtrix / Halftrix Accessories

OS-1 Hospital grade outlet strip. 15 foot long cord with 6 outlets.

CW-1 Cord wrap. Usually used with outlet strip (OS-1) sold separately.

UR-1, UR-2, UR-3 Utility receptacles, available in three sizes, UR-1 (0.85 liters), UR-2 (1.5 liters), UR-3 (2.8 liters). Requires accessory rail.

IV-2 Chrome IV pole and mounting bracket. Mounts on Unicart, Med Jr, Medtrix, DMC. For ER2000 please order part number IV-3.

FT-1 Fluid/equipment storage tray system for Steel, Med Jr or Medtrix/Halftrix carts. Includes three trays. Aluminum carts purchase FTA-1. Each tray measures 17.52"W x 4.04"D x 4.5"H

WNH-1 Wire notebook holder. Mounts to accessory railing.

GMH-1 Glove & Mask box holder. Made from lightweight aluminum. Adjusts to many sizes of boxes. Mounts to railing on cart.

CHT-1 Chart Holder. Mounts to accessory railing. 13"W x 4"D x 8.5"H

TRIM-MT

Custom color your cart to match you décor. We offer four colors of drawer trim to choose, Blue, Green, Pink, or Violet.

Please use the following part numbers. The last letter will specify color
Half-sized Drawers: TRIM-MTH(B)
Full-sized Drawers: TRIM-MTF(B)

B=blue, G=green, P=pink, V=Violet
TRIM-MTFB
TRIM-MTFG
TRIM-MTFP
TRIM-MTFV
TRIM-MTHB
TRIM-MTHG
TRIM-MTHP
TRIM-MTHV




WNH-1



GMH-1




CHT-1



TRIM-MT

Medtrx / Halftrx Cart Configurations

Set-up MED-TRX 1.5 (full size) and HALF-TRX 1.0 (compact) medication carts. The chart below indicates the required spaces for each item. (Side A is crossed out on the HALF-TRX 1.0 configuration).



Half size drawers.

Full size drawers and cassettes.

| Med-Trax Item # | Description | Spaces Required |
|-----------------|---------------------------|-----------------|
| DRW-3H | 3" Half Size Drawer | 2 |
| DRW-4H | 4" Half Size Drawer | 3 |
| DRW-7H | 7" Half Size Drawer | 5 |
| DRW-9H* | 9" Half Size Drawer | 6 2/3 |
| DRW-10H | 10" Half Size Drawer | 7 |
| DRW-3F | 3" Full Size Drawer | 2 |
| DRW-4F | 4" Full Size Drawer | 3 |
| DRW-6F | 6" Full Size Drawer | 4 |
| DRW-7F | 7" Full Size Drawer | 5 |
| DRW-9F* | 9" Full Size Drawer | 6 2/3 |
| DRW-10F | 10" Full Size Drawer | 7 |
| CASS-3MTC | 2 Tier Cassette w/3" bins | 4 |
| CASS-5MTC | 2 Tier Cassette w/5" bins | 4 |
| CASS-7MTC | 2 Tier Cassette w/7" bins | 4 |

Full Drawers: 16"D x 17.25"W
Half Drawers: 16"D x 9.125"W

*Must be ordered in sets of three.

MEDTRX 1.5

Dimensions (includes casters)
Overall shell size: 39"W x 20"D x 43"H

HALFTRX 1.0

Dimensions (includes casters)
Overall shell size: 27"W x 20"D x 43"H

A completed set-up sheet is required with orders on custom MED-TRX and HALF-TRX Carts. Please contact Customer Service or your Sales Representative if you require assistance with set-up. Custom carts are not returnable.

Storing narcotics?

Individual locks can be added to medication cart drawers (not available for 3"H). Drawers are keyed separately from cart to add extra security. Ask your sales representative or call our customer service department at 1-800-833-4419 for more information.



NL-MT Narcotics Lock

Terms & Conditions

ORDERING AND SHIPPING INFORMATION

WATERLOO HEALTHCARE • 3730 EAST SOUTHERN AVENUE • PHOENIX, AZ 85040
REMIT TO ADDRESS: P.O. BOX 53555 • PHOENIX, AZ 85072-3555

CALL US
8:00 – 5:00 CST

OFFICE

(800) 833-4419 FAX: (602) 437-2270

www.waterloohealthcare.com

sales@waterloohealthcare.com

PLACING YOUR ORDER

- When placing your order with us, please refer to your customer account number, which appears on all of our invoices and packing lists.
- Please write your customer account number on the first page of your catalog for immediate reference when placing orders
- All orders must be in writing. Please fax your order to us or send in by mail.
- Although they may not be listed individually, Customer Service can help you order special parts; please just ask.

TERMS OF SALE

- Net 30 days for regular customers with established credit.
- Visa, MasterCard, and American Express accepted.
- Service or handling charges are not refundable.

PRICING

- Pricing is subject to change without notice and products are priced "each" unless otherwise noted.
- When calling our Customer Service Department for any reason, please ask if your price listings are current.
- This catalog is available free of charge to any qualified individual, organization, or institution.

SHIPPING TERMS

- Our standard shipping terms are F.O.B. Phoenix, AZ.
- Stock items can be shipped within two days after receipt of your order.
- Wherever possible, we ship via UPS Ground Service. For emergency orders, we will gladly ship your order via UPS Air for guaranteed delivery within one or two days.
- All orders that require a pallet will ship via LTL Carrier

NATIONAL SALES FORCE

- Waterloo Healthcare supports our customers with a National Sales Force to assist with product demonstrations or any technical questions about our products. Please ask our Customer Service Department if there is a representative in your area.

RETURNED GOODS

- All returned goods must have an authorization number assigned by our Customer Service Department. Returns must be within 30 days from date of original shipment.
- Please refer to your purchase order number and customer account number when phoning in your request for returning merchandise.
- All returned goods must be sent prepaid. Any return that is not prepaid will be refused.
- Any cart with Electronic or Pushbutton locking systems, materials that have been modified or are not in resalable condition are non-returnable. Large or Custom orders are subject to contract and are non-returnable.

RESTOCKING CHARGE

- A 25% restocking charge may be charged to your account for returned merchandise.
- An exception to this policy would be if we shipped incorrect merchandise. If this is the case, please call our Customer Service Department and report the incorrect shipment. The same policy regarding an authorization number applies.

DAMAGED GOODS

- All damage must be reported to WHC the day the items are delivered.
- If you receive a shipment via truck, inspect all cartons at time of delivery. If you receive a package that you believe is damaged, we ask that you do not sign for it, open all packaging and if there is damage call 1-800-833-4419 as soon as possible. **YOU MUST INSPECT AND MAKE NOTE OF ALL possible damage on Carriers Delivery Receipt or Bill of Lading**
- You may be asked to take pictures of the damage and send them to WHC to file the claim.
- Filing claim for damage not marked on carriers delivery receipt at time of delivery (also know as "concealed damage") will be the responsibility of the customer.
- Please save all boxes and packing material to prove that the items were packed properly.

PRODUCT WARRANTIES

- Waterloo Healthcare (Seller) warrants the products it manufactures to be free from defects in material and workmanship under normal and proper use and service for a period not exceed five (5) years from the date of delivery to the original location.
- This warranty does not apply to any product that has been subject to abuse, misuse, negligence, modification, normal wear, or an accident (dents & scratches are considered normal wear).
- Electronic components and plastic parts have a one (1) year warranty under the aforementioned conditions of use.
- The Seller's liability is limited to the cost of the repair or replacement of any products (at its factory) which fail to comply with the foregoing warranty. In no event shall the seller be liable for any consequential damages claimed as a result of breach of the foregoing warranty.
- The Purchaser must reference the original purchase order number or Seller's invoice on any claims. The Seller will determine if the products should be returned to the factory or if parts should be sent to the customer for repair of the product.
- The warranty stated herein is in lieu of all warranties, expressed or implied, including but not limited to, merchantability or fitness for a particular purpose.

Index

Accessories

| | |
|------------------------|-------|
| Medication Cart | |
| MED-TRX | 52-56 |
| Dual Sided (DMC) | 52-56 |
| Unicart | |
| Aluminum | 32-39 |
| Bedside | 32-39 |
| Slim Med. Junior | 32-39 |
| Steel | 32-39 |

Accessory Packages

| | |
|-----------------------------|------------|
| Airway Breathing | 22, 40 |
| Anesthesia | 16, 17, 40 |
| Cast | 23, 40 |
| Emergency | 11, 12, 40 |
| Isolation/Infection | 18, 19, 40 |
| IV Therapy | 23, 40 |
| Pediatric | 15, 40 |
| Treatment & Procedure | 22-24, 40 |

Aluminum Carts

| | |
|---------------------------|----------------|
| Anesthesia | 16-17 |
| Dividers | 38-39 |
| Electronic Locking | 5-7, 16, 28-29 |
| Emergency | 11-15 |
| Features & Benefits | 6-7 |
| Isolation/Infection | 18 |
| Locking Systems | 6 |
| Medication | 43-45 |
| Pediatric | 15 |
| Proximity Locking | 5-6 |

Ambu-Bag Hooks

Anesthesia Unicart

| | |
|---------------------------------|----------------|
| Aluminum | 16 |
| Electronic Locking | 16 |
| Steel | 17 |
| Automated Locking Systems | 4-7, 16, 28-29 |

Backboard *See Cardiac Board & Bracket*

Bar-Accessories *See Side Rail or Cross Bar*

Base Rail System

| | |
|----------------------|--------|
| Aluminum Carts | 16, 35 |
| Med. Junior | 35 |
| Steel Carts | 17, 35 |

Basket

| | |
|---------------------------|----|
| Utility Receptacles | 34 |
| Waste | 33 |
| Wire | 34 |

Bins

| | |
|--------------------|------------------|
| Tilt Bins | 34 |
| Patient Bins | 43, 45-49, 52-53 |

Breakaway Locking System

| | |
|---------------------------|-----------------------|
| Aluminum Unicarts | 6-7, 11-15, 18, 22-24 |
| Bedside | 4, 20-21 |
| Emergency Carts | 10, 11-15 |
| Gatlock | 37 |
| Med Junior Unicarts | 4, 19, 20-21 |
| Steel Unicarts | 8-9, 11-15, 19-23 |

| | |
|--------------------------------|----|
| Breakaway seals (PS-100) | 37 |
|--------------------------------|----|

Cabinets

| | |
|-----------------|----|
| Narcotics | 42 |
|-----------------|----|

Cardiac Board & Bracket

| | |
|-------------------------|----|
| Aluminum Unicarts | 32 |
| Economical | 32 |
| ER2000 | 32 |
| Steel Unicarts | 32 |

Cart

| | |
|-------------------------------|-------------------------------------|
| Airway | 22 |
| Aluminum Unicarts | 5-7, 11-16, 18, 23-24, 28-29, 43-45 |
| Anesthesia | 16-17 |
| Bedside | 4, 20-21 |
| Med Junior | 4, 14, 19-21 |
| Medication | 41-51 |
| Mobile Workstation | 24 |
| Simulated Wood Cart | 24, 44 |
| Steel Unicarts | 8-9, 11-15, 17, 19-27 |
| Split Drawer Cart | 22 |
| Treatment and Procedure | 22-27 |

Cassette Mounting Bracket

Cassettes (Medication Carts)

| | |
|--------------------|---------------|
| Dual Sided | 43, 47, 49 |
| Single Sided | 43, 46-48, 52 |

Cast Cart

Catheter Holder

Chart Holder

Coat Rack

Code Cart *See Emergency or ER2000*

Color Options

| | |
|----------------------------|---------------------------------|
| Med Junior | 4, 20-21 |
| Medication Cart Trim | 56 |
| Pediatric | 15 |
| Unicarts | 4-5, 7, 9, 13, 15, 18-20, 22-24 |

Cord Wrap

CPR Backboard *See Cardiac Board & Bracket*

Crash Cart *See Emergency*

Cross Bar

Cup Dispenser

Defibrillator Shelf *See Swivel Monitor Shelf*

Dividers

| | |
|-------------------------|-------|
| Aluminum Unicarts | 38-39 |
| Bedside | 39 |
| Med. Junior | 39 |
| Steel Unicarts | 38-39 |

Medication

| | |
|--------------------|-------|
| Dual Sided | 53 |
| Single Sided | 53-54 |

Dual Sided Medication Carts

| | |
|---------------------------|-------|
| Features & Benefits | 49 |
| Patient Bin Carts | 47-49 |

Economy Carts

Electronic Locking Systems

| | |
|---------------------------|----------------|
| Aluminum Unicarts | 5-7, 16, 28-29 |
| Med Junior Unicarts | 20-21 |

Index

Emergency Carts

| | |
|-----------------|-------|
| Aluminum | 11-15 |
| ER2000 | 10 |
| Pediatric | 15 |
| Steel | 11-15 |

| | |
|---------------------|----|
| ER2000 | 10 |
|---------------------|----|

Exchange Cassettes *See Cassettes*

| | |
|---------------------------------|----|
| Gatelocking System | 37 |
|---------------------------------|----|

| | |
|---------------------------|----|
| Glove Holder | 33 |
|---------------------------|----|

HALFTRX Carts

| | |
|---------------------------|--------|
| Dividers | 53, 54 |
| Features & Benefits | 46 |
| Locking Systems | 46 |
| Patient Bin Carts | 48 |
| Pill Carts | 51 |

Hooks

| | |
|----------------------|----|
| Utility Hooks | 32 |
| Ambu-Bag Hooks | 32 |

Infant Color Coded Cart *See Pediatric*

Instant Access Emergency Cart *See ER2000*

Isolation/Infection Unicart

| | |
|-----------------------|----|
| Aluminum | 18 |
| Slim Med Junior | 19 |
| Steel | 19 |

| | |
|----------------------|----|
| IV Pole | 33 |
|----------------------|----|

Locking Security Options

| | |
|-------------------------|--------------|
| Aluminum Unicarts | 5-6, 37 |
| Gatelock | 37 |
| Medication Carts | 44-46, 54-55 |
| Steel Unicarts | 8, 37 |

Labels

| | |
|-----------------------|----|
| ER2000 | 10 |
| Medication Bins | 53 |

Laptop Shelf *See Swivel Monitor Shelf*

Lever Locking Systems

| | |
|-------------------------|-----------------------|
| Aluminum Unicarts | 6-7, 11-15, 18, 23-24 |
| Steel Unicarts | 8-9, 11-15, 19-24 |

Lightweight Cart *See Aluminum*

Locking Systems

| | |
|-------------------------|--------|
| Aluminum Unicarts | 6 |
| Gatelock | 37 |
| Medication Carts | 44, 46 |
| Steel Unicarts | 8 |

| | |
|--------------------------|----|
| Mask Holder | 33 |
|--------------------------|----|

Med Junior

| | |
|---------------------------|----------|
| Bedside | 4, 20-21 |
| Dividers | 39 |
| Emergency | 14 |
| Features & Benefits | 4, 20-21 |
| Isolation/Infection | 19 |
| Proximity Locking | 4, 20-21 |

Medication Cart

| | |
|---------------------------|--------|
| Dividers | 52-54 |
| Dual Sided (DMC) | |
| Features & Benefits | 47, 49 |
| Locking Systems | 47, 49 |

Single Sided

| | |
|----------------|----------------|
| Aluminum | 43-45 |
| HALFTRX | 43, 46, 48, 51 |
| MEDTRX | 43, 46, 48, 51 |

MEDTRX Carts

| | |
|---------------------------|--------|
| Dividers | 52, 54 |
| Features & Benefits | 46 |
| Locking Systems | 46 |
| Patient Bin Carts | 48 |
| Pill Carts | 51 |

Narcotics

| | |
|----------------------|--------------------|
| Cabinet | 42 |
| Gatelock | 37 |
| Locking drawer | 43, 46, 48, 51, 54 |

Needle Container *See Sharps Bracket*

Notebook Holder *See Chart or Wire Notebook*

| | |
|--|-------|
| Organizers, Clear Plastic | 38-39 |
|--|-------|

| | |
|---------------------------|----|
| Outlet Strip | 36 |
|---------------------------|----|

| | |
|-----------------------------|----|
| Overhead Shelf | 35 |
|-----------------------------|----|

| | |
|----------------------------------|----|
| Oxygen Tank Bracket | 32 |
|----------------------------------|----|

Packages *See Accessory Packages*

Patient Bins (Medication Carts)

| | |
|--------------------|----------------|
| Dual Sided | 43, 47, 49, 53 |
| Single Sided | 43, 46, 48, 53 |

Pediatric Unicart

| | |
|----------------|----|
| Aluminum | 15 |
| Steel | 15 |

| | |
|------------------------------|----|
| Phlebotomy Cart | 25 |
|------------------------------|----|

Pill Card Carts *See Medication Carts*

| | |
|---|----|
| Pill Crusher & Pouches | 55 |
|---|----|

| | |
|----------------------------|----|
| Plastic Seals | 37 |
|----------------------------|----|

| | |
|-------------------------------|----|
| Portable Storage | 41 |
|-------------------------------|----|

Pouches *See Pill Crusher*

Power Strip *See Outlet Strip*

Procedure Cart *See Treatment & Procedure*

Proximity Locking Carts

| | |
|------------------|------------------------------|
| Aluminum | 5-7, 16, 23-24, 28-29, 43-45 |
| Bedside | 4, 20-21 |
| Med Junior | 4, 20-21 |

Push Button Locking Systems

| | |
|---------------------------|----------|
| Aluminum Unicarts | 6-7 |
| Bedside Unicarts | 4, 20-21 |
| Med Junior Unicarts | 4, 20-21 |
| Medication Carts | 46 |
| Steel Unicarts | 8-9, 17 |

Railing & Supports

| | |
|-------------------|----|
| Aluminum | 35 |
| Med. Junior | 35 |
| Steel | 35 |

| | |
|--|----|
| Seals, Plastic breakaway (PS-100) | 37 |
|--|----|

Security Cabinet *See Narcotics Cabinet*

Security Drawer

| | |
|------------------------|----|
| Aluminum | 37 |
| Medication Carts | 54 |
| Steel Unicart | 37 |

| | |
|--|----|
| Sharps Bracket (Adjustable) | 32 |
|--|----|

Index

| | |
|--|----------|
| Sharps Bracket (Lockable) | 32 |
| Shelf(s) | |
| Cast | 23, 33 |
| Overhead | 35 |
| Side | 35 |
| Suction | 35 |
| Suture | 35 |
| Swivel | 35 |
| Side Storage (DMC)..... | 47 |
| Side Rail (Steel Unicart accessory Bar) | 36 |
| Slim “Med. Junior” Cart | |
| Bedside | 4 |
| Dividers | 39 |
| Emergency | 14 |
| Features & Benefits | 4, 20-21 |
| Isolation/Infection | 19 |
| Proximity Locking | 4, 20-21 |
| Specimen Collection <i>See Phlebotomy</i> | |
| Steel Unicart | |
| Anesthesia | 17 |
| Bedside | 4, 20-21 |
| Cast | 23 |
| Dividers | 38-39 |
| Emergency | 11-15 |
| Features & Benefits | 8-9 |
| Isolation/Infection | 19 |
| Locking Systems | 8-9 |
| Pediatric | 15 |
| Phlebotomy | 25 |
| Med Junior | 4, 20-21 |
| Slim | 4, 20-21 |
| Treatment & Procedure | 22-27 |
| Storage Basket | |
| Tilt Bins | 34 |
| Utility Receptacles | 34 |
| Waste | 33 |
| Wire | 34 |
| Suction Shelf | 35 |
| Suture Shelf | 35 |
| Swivel Monitor Shelf | 35 |
| Tape Dispenser | 34 |
| Terms & Conditions | 58 |
| Tilt Bins | 34 |
| Top, replacement (Unicart) | 36 |
| Trash Receptacle <i>See Waste Basket</i> | |
| Tray, Removable <i>See Divider</i> | |

| | |
|---|-------------------------------------|
| Treatment & Procedure Unicart | |
| Airway | 22 |
| Aluminum | 22-24 |
| Cast Cart | 23 |
| IV Therapy | 23 |
| Mobile Workstation | 24 |
| Phlebotomy | 25 |
| Simulated Wood Cart | 24 |
| Trellis System <i>See Base Rail System</i> | |
| Trolley <i>See Carts</i> | |
| Unicart | |
| Aluminum | 5-7, 11-16, 18, 23-24, 28-29, 43-45 |
| Economy | 26-27 |
| Med Junior | 4, 19-21, 26, 30 |
| Steel | 8-9, 11-15, 19-27, 28, 30 |
| Utility | |
| Hooks | 32 |
| Receptacles..... | 33 |
| Wall Lock | 55 |
| Waste Basket | 33, 55 |
| Wire | |
| Accessory Basket | 34 |
| Notebook Holder | 34 |
| Storage Basket | 34 |

Color Options



GRA-Charcoal



DKB-Dark Blue



ELB-Electric Blue



FWG-Fairway Green



LTB-Light Blue



LTG-Light Gray



NTL-Natural



ORG-Orange



PNK-Pink



RED-Red



SPG-Spring Green



TLG-Teal



VIL-Violet



WHT-White



WNE-Wine



YEL-Yellow

Cart Specifications

| Overall Dimensions | Exterior Colors See page 62 | Drawer Colors See page 62 | Drawer Sizes | Lock Options | Weight Range (lbs) | Built-in Dividers | Antimicrobial | Pull-Out Shelf |
|--------------------------------------|--------------------------------|------------------------------|------------------|---|-----------------------|----------------------|---------------|-------------------|
| STEEL UNICART | | | | | | | | |
| Short 37.5"H x 29"W x 24.5"D | LTG,RED,ELB | All 16 options | 3,4,6,9,12,Panel | Lever, Key, Pushbutton | 85-100 | No | Yes | No |
| Tall 43.5"H x 29"W x 24.5"D | LTG,RED,ELB | All 16 options | 3,4,6,9,12,Panel | Lever, Key, Pushbutton | 90-120 | No | Yes | No |
| X-Tall 46.5"H x 29"W x 24.5"D | LTG,RED,ELB | All 16 options | 3,4,6,9,12,Panel | Lever, Key, Pushbutton | 100-140 | No | Yes | No |
| ALUMINUM UNICART | | | | | | | | |
| Short 36"H x 29"W x 24.5"D | LTG,RED,ELB | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 65-90 | Yes | Yes | Yes |
| Medium 39"H x 29"W x 24.5"D | LTG,RED,ELB | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 70-95 | Yes | Yes | Yes |
| Tall 42"H x 29"W x 24.5"D | LTG,RED,ELB | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 75-95 | Yes | Yes | Yes |
| X-Tall 45"H x 29"W x 24.5"D | LTG,RED,ELB | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 80-100 | Yes | Yes | Yes |
| STEEL BEDSIDE / MED JR | | | | | | | | |
| Short 36"H x 23"W x 24.5"D | Light Gray | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 80-95 | No | Yes | No |
| Tall 42"H x 23"W x 24.5"D | Light Gray | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 85-115 | No | Yes | No |
| X-Tall 45"H x 23"W x 24.5"D | Light Gray | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 90-130 | No | Yes | No |
| ALUMINUM BEDSIDE / MED JR | | | | | | | | |
| Short 36"H x 23"W x 24.5"D | Light Gray | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 60-80 | Yes | Yes | Yes |
| Medium 39"H x 23"W x 24.5"D | Light Gray | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 60-85 | Yes | Yes | Yes |
| Tall 42"H x 23"W x 24.5"D | Light Gray | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 65-85 | Yes | Yes | Yes |
| X-Tall 45"H x 23"W x 24.5"D | Light Gray | All 16 options | 3,6,9,12 | Lever, Key, Pushbutton, Electronic, Proximity Card | 70-90 | Yes | Yes | Yes |

Need a custom cart that doesn't meet any of the above? Just call us!

Cart Specs



*Waterloo
Healthcare*

Waterloo Healthcare
3730 E. Southern Avenue
Phoenix, AZ 85040-3958

CHANGE SERVICE REQUESTED