



Perfecseal, a Bemis Company, 3500 N Main St, Oshkosh, WI, 54901
www.perfecseal.com 920-303-7000 800-633-7587

Medi-Oxide® Tyvek® Self Seal Pouches



Medi-Plus® brand sterilization packaging is manufactured by Perfecseal the largest manufacturer of medical grade pouches in the world. Our Oshkosh, WI manufacturing facility is ISO 13485 and FDA registered. The quality systems and good manufacturing processes employed under these certifications assure your institution of receiving products equal in construction and quality to the products supplied to our medical device manufacturer customers.

Medi-Oxide sterilization pouches and tubing are indicated for low temperature sterilization methods such as hydrogen peroxide gas plasma, ethylene oxide gas, and peracetic acid. All films used in our pouches and tubing are engineered specifically for medical applications. These medical grade films sealed to paper or Tyvek® provide a unique combination of superior seal strength, clean consistent peel, and sterility assurance.

Medi-Oxide® Tyvek® self-seal pouches feature the original pre-fold – assuring accurate, consistent and fast sealing without perforation holes. Unique angled graphics ease differentiation from paper pouches to prevent accidental use in high temperature steam autoclaves.

Call (800) 633-7587 for your FREE trial box.

Reorder Number	Size (inches)	Size (metric)	No. Inner	Units Inner	Total
98-54826	3 ½" x 8 ½"	90mm x 216mm	10	100	1000
98-54828	5 ¼" x 10 ¼"	133mm x 260mm	10	100	1000
98-54832	5" x 15"	127mm x 381mm	10	100	1000
98-54829	7 ½" x 13 ½"	191mm x 343mm	10	100	1000
98-54830	12" x 15 ½"	305mm x 394mm	5	100	500
98-54834	4" x 24"	102mm x 610mm	10	100	1000

For Healthcare Professionals Who Care from Healthcare Packaging Professionals Who Care

Tyvek® is a registered trademark of the DuPont Company