



Technical Information
 General Procedure Cart
 Surgical, Procedure, Treatment

SLIM MED JR SERIES CART
4 DRAWER

Model Shown:
 JSGKU-3369-LTG

Designed with input from Healthcare Professionals this cart provides controlled access to critical supplies, narcotics and other medications which must be secured to ensure availability in an urgent situation.

Built to last this with a five (5) year bumper to bumper limited warranty (see terms and conditions). This Short, Light Gray / Light Gray, 4-drawer cart comes fully assembled and includes:

- o Antimicrobial properties built in to protect majority of the product.
- o Removable colored drawer fronts
- o Pre-threaded holes for accessory placement
- o Dual ergonomic push handle
- o Drawers with ball bearing full extension slides
- o Centralized locking with 5 available options! (Breakaway, Key, Pushbutton, Electronic, Proximity)
- o 5" premium casters (2 braking, 1 tracking, 1 swivel)
- o ABS plastic molded top.
- o Stabilizer base and high impact plastic bumper



Dimensions:


- Product: 23"W x 24.5"D x 36"H
- Drawers: (2) 3"H, (1) 6"H, (1) 9"H
- Shipping: 32"W x 32"D x 45"H

Approximate Weights:

- Product: 84 lbs
- Shipping: 95 lbs


Built to order every time! **CUSTOM CONFIGURATIONS** available on several styles **with NO extra fee!**

Exterior Color Options: **Antimicrobial properties**



LTG-Light Gray

Drawer Front Color Options:



Customizable drawer configurations!



View our full online catalog at www.waterloohealthcare.com
 Contact customer service at: 1-800-833-4419
 Email: Sales@waterloohealthcare.com

This information is subject to change without formal notice

Waterloo Healthcare, Inc.