

Medication Cart Systems

Waterloo Healthcare medication carts (trolleys) offer quality...security...and...flexibility.



MEDTRX 1.5

Dimensions (includes casters)
Overall exterior size 39"W x 20"D x 43"H



MEDTRX-3AC

*Shown with optional narcotics drawer lock (NL-MT)
Pictured as part #MTK14B771/TFTFT*

MEDTRX 1.0

Dimensions (includes casters)
Overall exterior size 27"W x 20"D x 43"H



HALFTRX-2LTC

*Pictured as part #HTK1-999
Holds 450 pill cards.*

CASSETTE UNICART



CASSETTE UNICART

*Pictured above as part
#UTGKU-C3336-LTG*

DMC-3 (includes side module)

Dimensions (Includes stabilizer base, side handle and casters)
Overall32"W x 24.5"D x 43.5"H



DMC3-4S10

PHOENIX LTC CARTS

Overall exterior size: 38"W x 24"D x 43"H



750PC
Holds 750 pill cards.

Medication Cart

The NEW Aluminum Wood Simulated Cart from Waterloo Healthcare



Features & Benefits

- Antimicrobial properties built-in to protect portions of the product
- Elite Aluminum offers easier maneuverability and safety with a lightweight aluminum construction being 30% lighter than standard steel carts
- ABS plastic molded top.
- New dual pull-out shelves provide additional work surface on both sides of the cart.
- Adjustable dual push handles.
- Centralized locking with 5 available options.
- Drawers with ball bearing slides designed to close and lock even after the cart has been secured.
- Drawers include built-in, interlocking dividers corresponding to the drawer height.
- 5" Premium easy roll swivel casters (2 brake, 1 tracking, 1 full swivel).
- Stabilizer base with high impact plastic bumper to add stability.
- Latex Free.
- Available in one size:
Tall: (27" drawer height)
29" W x 24.5" D x 42" H
and 12-inch deep drawers.



Simulated Wood Cart

For a less industrial feel choose the Simulated Wood cart to match the décor of your facility.

Customize to use for unit dose boxes or use as a pill card cart with the ability to configure with 3 inch, 6 inch, 9 inch or 12 inch deep drawers.

UTWKA-33669-WOD

Lightweight Aluminum Unicart with wood veneer.

Overall Dimensions: 29"W x 24.5"D x 42"H
(Include casters, stabilizer base and ABS top)
Drawer Dimensions: 22"W x 16.5"D
Drawer Height: (2) 3"H, (2) 6"H, (1) 9"H

Weight: 74 lbs.

Aluminum Locking Options
Page 6

Aluminum Accessory Options
Pages 32-37

Self Configuration - Tall Only
Page 29

UNIT DOSE CART

UXGSA-333333336-LTG

(Cart shown with optional pushbutton locking system).

Dimensions

(includes casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 45"H

Drawers: 22"W x 16.5"D

Drawer Height: (8) 3"H, (1) 6"H

Weight: 85 lbs

5 available locking systems. See Page 6.



UXGSA-333333336-LTG



UXGKA-3999-LTG

PILL CARD CART

UXGKA-3999-LTG

Dimensions

(includes casters, stabilizer base and top)

Overall: 29"W x 24.5"D x 45"H

Drawers: 22"W x 16.5"D

Drawer Height: (1) 3"H, (3) 9"H

Weight: 85 lbs

5 available locking systems. See Page 6.

Holds 450 pill cards.



DIV-UC9

DIV-UC9 Divide the drawer into rows to provide organization for pill cards or fluids. For additional Aluminum cart accessories see page 32 - 39.



NARC-BOX

NARC-BOX Secure medications in the 9"H aluminum or steel Unicart drawers.

Larger View – See Page 62



GRA-Charcoal



DKB-Dark Blue



ELB-Electric Blue



FWG-Fairway Green



LTB-Light Blue



LTG-Light Gray



NTL-Natural



ORG-Orange



PNK-Pink



RED-Red




SPG-Spring Green



TLG-Teal



VIL-Violet



WHT-White



WNE-Wine



YEL-Yellow

Medication Cart Highlights

Waterloo Healthcare medication carts offer quality, security and flexibility. With two locking options and flexible custom configuration, the Medtrix-series carts will fit the needs of any facility. See page 57 for customer configuration options.

Our Medtrix series is available in two sizes. Our full size cart offers half size drawers on the left for storage and supplies and space for up to 50 patient bins or full size drawers on the right. We also have a half size cart (Halftrix) for use where space is limited. This cart can hold full size storage drawers, patient bins or a mix of both.

Locking options to choose from:



Tamper Resistant Best® Key Locking System



Simplex Pushbutton With Key Lock Override

PATIENT CASSETTE CART SYSTEMS

Styled push handles

Tamper resistant Best® Lock

Open/closed warning label on both sides

Optional narcotic lock (not available on 3"H drawers)



MEDTRX-1AC

Pictured as part #MTK1-4B411/6TTTT

PILL CARD CART SYSTEMS

Antimicrobial properties built-in to protect portions of the product.



Antimicrobial Surfaces



HALFTRX-1LTC

Pictured as part #HTK1-33311

Slide out shelf

ABS plastic removable molded top

Pill card drawers hold a total of 300 pill cards.

Wrap around bumper

5" premium easy roll casters (one swivel, one tracking, two w/locking toe brake)

*Quantity will vary depending on card thickness and use of dividers.

Medication Cart Highlights

When dual-sided accessibility is required, look no further, a time proven design from Waterloo Healthcare, our dual-sided medication carts have been popular worldwide for many years.

PATIENT CASSETTE CART SYSTEMS

Pull-out shelf DMC3-4S10



Both sides of cart have same drawer configuration.

- Cup holder
- Push handle
- Positive center catch – prevents drawers from accidentally passing through
- Lockable/portable cassette storage
- Multiple drawer sizes
- Stabilizer Base with Bumper Guard adds stability by lowering the center of gravity on the cart, keeping it from tipping
- 5” premium easy roll casters (one swivel, one tracking, two w/locking toe brake)

Antimicrobial properties built-in to protect portions of the product.




DMC-3 Dual-Sided Medication Cart
 Dimensions (includes casters)
 Overall32”W x 24.5”D x 43.5”H

SIDE VIEW



DRF-2

Add lockable security to storage area on side module. Ideal for liquid storage.



- Removable, replaceable ABS top
- Tamper resistant Best® Lock
- Storage shelves
- Dual access drawers
- Roller bearing slides
- Tilt waste receptacle

Medication Cart