



Change Your Vue . . .

for laparoscopic & robotic-assisted laparoscopic surgeries



replace all this . . .

with LaparoVue™ Visibility System
an All-in-One Solution



Pending 510K approval.

Product Ordering Information

Description	Catalog Number	Quantity
LaparoVue™ Visibility System	LAPVUE10	10/Case

To view additional **Laparoscopic Plume Filtration Devices** and other solutions visit www.buffalofilter.com.

LaparoVue™ Visibility System



All-in-One Solution
for improved laparoscopic visibility

Buffalo Filter, LLC
5900 Genesee Street, Lancaster, NY 14086 U.S.A.
1.716.835.7000 1.800.343.2324

PN 904424 REV A Actual product may differ slightly in appearance to images shown.
* Photographs are representative of a fogged and clear / defogged lens but do not reflect actual product use.
Buffalo Filter® is a registered trademark of Buffalo Filter, LLC in the United States and other countries. © 2015 Buffalo Filter, LLC. All Rights Reserved.



Laparoscopic Visualization ...

Complexity Made Simple



warms & white balances

- warms scope in 5 minutes to prevent fogging
- unit warms up to 5 hours
- allows for true white balance




cleans

- provides crisp, clear view in 5 seconds
- minimizes procedural delays
- radiopaque VueTip™ Trocar Swabs remove debris and excess liquid from trocar cannula
- x-ray detectable Microfiber Cloth designed to capture residue and particulates at a microscopic level

defogs

- pre-filled cleaning / defogging solution
- helps prevent intra-operative scope fogging



Distinctive By Design ... *unique product attributes*

one unit, two ports

- two ports: warming / white balancing and cleaning / defogging
- fits 3mm–12mm scopes without port size converter
- allows for simultaneous use of both ports

detachable scope cradle

- safely and securely supports scope
- reduces tension on lens and potential for damage

built-in LED lights

- provides enhanced port visibility

ergonomic design

- stable base prevents tipping
- contoured, easy grip profile
- engineered with soft silicone material

